

Philips Avent  
Twin electronic breast  
pump



SCF294/02

# Designed for Comfort

Inspired by nature

Our unique pump is designed with an electronic memory that puts you in control by learning and continuing your personal pumping style.

#### **Benefits and features:**

- Ultra comfortable and effective
- Control at your fingertips
- Infinitely variable speed, vacuum and rhythm – no pre-sets
- Quiet, effortless and portable
- Pumps more milk faster

#### **Stimulate and Express**

- Soft massage cushions trigger natural let-down, like baby
- Gentle vacuum mimics baby's suckling for steady milk flow

**PHILIPS**  
**AVENT**

# Highlights

## Soft massage cushions



Patented massage cushions of the Philips Avent breast pump flex in and out, replicating baby's suckling action and is intended to help stimulate a fast natural let down.

## Gentle-draw vacuum

The gentle-draw vacuum mimics baby's suckling to create a steady milk flow that requires less pumping.

## Patented Let-down Massage Cushion



Five, soft petal massagers gently flex in and out working in tandem with the reliable vacuum. This is intended to help stimulate a fast, natural milk flow.

## Pumps more milk faster



Double pumping has been shown to improve milk production. The twin electronic breast pump is the only breast pump designed to maintain the same level of vacuum and control for both pumps, which helps to ensure effective use.

## Infinitely variable control

Enjoy complete control over speed, vacuum and pattern – you can even set the pump to match your baby's own suckling action.

## Intuitive control handle



Everything that controls the pump is at your fingertips and close to the breast, which means you can make adjustments instantly and intuitively- no need to reach elsewhere.

## Quiet, effortless and portable



Quiet, effortless and portable

# Specifications

## Country of origin

- England

## Development stages

- Stage: 0 - 6 months

## Power

- Voltage: 100 - 240 V

## What is included

- Breast milk container (125 ml/ 4 oz): 2 pcs
- Extra soft newborn flow nipple: 2 pcs
- Bottle stand/funnel cover: 4 pcs
- Sealing disc for milk storage: 2 pcs
- Newborn nipple travel pack: 2 pcs
- Cool pack: 4 pcs
- Spare parts: 1 pcs
- Breast milk container (260 ml/ 9 oz): 2 pcs

## What is included

- Twin Electronic Breast pump: 1 pcs
- Manual pump parts for when power is not available: 1
- Microfiber travel bag: 1 pcs
- Thinsulate™ carrier: 2 pcs



Issue date 2021-02-15

Version: 7.3.1

© 2021 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)