

# PHILIPS

Infrared lamp



HP3614



## Relieves muscular pains and common colds

Deep, focused warmth

Philips infrared lamp gives pain relief for muscles and joints. Its 100 Watt infrared warmth is comfortable and penetrates deep into skin, improving blood circulation thereby giving relief for focused areas of 20x30 cm.

**Philips - the light expert**

- More than 15 million infrared lamps sold

# Highlights

## Philips - the light expert

Philips is one of the global leading companies in light, with over 100 years of light experience. Globally Philips already sold more than 15 million of Infrared lamps.

# Specifications

## Easy to use

Adjustability: 0-40 degrees (backwards)

## Weight and dimensions

Product dimensions: 20.0x12.0x16.0

(HxWxD) cm

Product weight: 0.5 kg

## Technical specifications

Power: 100 W

Voltage: 220/230 V

Cord length: 1.6 m

Frequency: 50 Hz

## Weight and dimensions

A-box dimensions: 22.4x26.2x58.8

(HxWxD) cm

F-box weight: 1 kg

Qty. on Euro pallet: 288 pcs

No. of F-boxes in A-box: 6 pcs

A-box weight: 3.1 kg

F-box dimensions: 21.1x12.6x19.2 (HxWxD) cm

## Technical specifications

Insulation: Class II (double isolation)

Type of lamps: R95, 100W

Lifetime of lamp: 750 sessions of 10 minutes

## Safety

IEC certified: In line with IEC 60335

## Logistic data

Country of origin: Germany

## Logistic data

CTV code: 884361401000

## Easy storage

Cord storage

