

# PHILIPS

 Product Picture

## AVENT

Single electronic  
breast pump

Out and about set

SCF292/13

## Designed for Comfort

Inspired by nature

Our unique Out & About BPA-Free pump is designed with an electronic memory that puts you in control by learning and continuing your personal pumping style.

### **Benefits and features:**

- Learns and continues your pumping rhythm
- With insulated travel bag and cool packs

### **Helps you breastfeed longer**

- More milk, naturally

### **Stimulate and Express**

- Soft massage cushions trigger natural let-down, like baby
- Gentle vacuum mimics baby's suckling for steady milk flow

### **Use three ways**

- By mains power, battery power and manual

# Highlights

## With insulated travel bag and cool packs



Everything you need to maintain your milk supply and to express and store your breast milk when away from your baby.

## Soft massage cushions



Patented massage cushions of the Philips Avent breast pump flex in and out, replicating baby's suckling action and is intended to help stimulate a fast natural let down.

## Gentle-draw vacuum

The gentle-draw vacuum mimics baby's suckling to create a steady milk flow that requires less pumping.

## Learns and continues your pumping rhythm

Electronic memory learns and at the touch of a button, continues your personal pumping rhythm.

## Mains, battery and manual

In addition to mains power, this single electronic pump can also operate using battery power and as a manual pump, for total flexibility.

## More milk, naturally



Expression allows breast milk to be available, even when you can't be there

# Specifications

## Country of origin

England

## Development stages

Stage: 0 - 6 months

## Power

Voltage: 100-240 V

## What is included

Bottle stand/funnel cover: 2 pcs

Breast milk container (125 ml/ 4 oz): 2 pcs

Newborn nipple travel pack: 1 pcs

Cool pack: 2 pcs

Spare parts: 1 pcs

Breast milk container (260 ml/ 9 oz): 1 pcs

## What is included

Battery pack: 1 pcs

Manual pump parts for when power is not available: 1

Microfiber travel bag: 1 pcs

Power cord: 1 pcs

