



Philips  
InfraCare infrared lamp

**200W**  
Focused treatment area



HP3621

## Effective pain relief

Deep penetrating infrared warmth

Deep penetrating warmth from infrared lamp

### Effective pain relief

- Deep penetrating Infrared warmth

# PHILIPS

# Specifications

## Easy to use

- Adjustability: 0-40 degrees (backwards)
- On/off switch
- Soft-touch handgrip: for easy transport and positioning

## Weight and dimensions

- Product dimensions: 20x30.9x17.2 (WxHxD) cm
- Product weight: 1.2 kg

## Technical specifications

- Power: 200 W
- Voltage: 110 or 220-230 V
- Cord length: 200 m
- Frequency: 50/60 Hz

## Weight and dimensions

- A-box dimensions: 29x40.4x74.5 (WxHxD) cm
- F-box weight: 1.6 kg
- A-box weight: 4.8 kg
- No. of F-boxes in A-box: 3 pcs
- Qnt. on Euro pallet: 60 pcs
- F-box dimensions: 28x38.6x24.5 (WxHxD) cm

## Focused treatment

- Treatment area: 20x30 cm
- Infrared halogen lamp: 200 W

## Technical specifications

- Insulation: Class II (double isolation)
- Lifetime of lamps: 500 hour(s)

## Safety

- IEC certified: Complies to IEC 60601 and 60335

## Logistic data

- Country of origin: Germany

## Medical appliance

- Medical Device Directive: 2007/47/EC, MDD 93/42/EEC

## Logistic data

- CTV code: 884362101000

## Effective pain relief

- Deep penetrating warmth

# Highlights

## Deep penetrating Infrared warmth



Infrared light has been proven to effectively relieve aches and pains caused by muscular pains and stiff joints. The comfortable warmth created by infrared light penetrates deep into your skin, stimulates blood circulation and warms your muscles. As your muscles are soothed, they automatically loosen up and relax. Because it makes the tissue more pliable, warmth also reduces stiffness and makes joints more flexible.



Issue date 2021-02-15

Version: 8.1.1

© 2021 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)