




Philips
InfraCare infrared lamp

200W
Focused treatment area



HP3621

Effective pain relief

Deep penetrating infrared warmth

Deep penetrating warmth from infrared lamp

Effective pain relief

- Deep penetrating Infrared warmth

PHILIPS

Specifications

Easy to use

- Adjustability: 0-40 degrees (backwards)
- On/off switch
- Soft-touch handgrip: for easy transport and positioning

Focused treatment

- Infrared halogen lamp: 200 W
- Treatment area: 20x30 cm

Safety

- IEC certified: Complies to IEC 60601 and 60335

Effective pain relief

- Deep penetrating warmth

Technical specifications

- Voltage: 110 or 220-230 V
- Power: 200 W
- Cord length: 200 m
- Frequency: 50/60 Hz

- Insulation: Class II (double isolation)
- Lifetime of lamps: 500 hour(s)

Weight and dimensions

- Product dimensions: 20x30.9x17.2 (WxHxD) cm
- A-box dimensions: 29x40.4x74.5 (WxHxD) cm
- Product weight: 1.2 kg
- F-box weight: 1.6 kg
- A-box weight: 4.8 kg
- Qnt. on Euro pallet: 60 pcs
- No. of F-boxes in A-box: 3 pcs
- F-box dimensions: 28x38.6x24.5 (WxHxD) cm

Logistic data

- Country of origin: Germany
- CTV code: 884362101000

Medical appliance

- Medical Device Directive: 2007/47/EC, MDD 93/42/EEC

Highlights

Deep penetrating Infrared warmth



Infrared light has been proven to effectively relieve aches and pains caused by muscular pains and stiff joints. The comfortable warmth created by infrared light penetrates deep into your skin, stimulates blood circulation and warms your muscles. As your muscles are soothed, they automatically loosen up and relax. Because it makes the tissue more pliable, warmth also reduces stiffness and makes joints more flexible.



Issue date 2022-07-28

Version: 8.1.1

EAN: 08 71010 32457 42

© 2022 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com