

# PHILIPS

 Product Picture

## AVENT

Twin electronic  
breast pump

SCF304/02

## Designed for comfort

Inspired by nature

Our unique BPA-Free pump is designed with an electronic memory that puts you in control by learning and continuing your personal pumping style.

### **Benefits and features:**

- Ultra comfortable and effective
- Control at your fingertips
- Infinitely variable speed, vacuum and rhythm – no pre-sets
- Quiet, effortless and portable
- Pumps more milk faster

### **Stimulate and Express**

- Soft massage cushions trigger natural let-down, like baby
- Gentle vacuum mimics baby's suckling for steady milk flow

# Highlights

## Gentle-draw vacuum

The gentle-draw vacuum mimics baby's suckling to create a steady milk flow that requires less pumping.

## Soft massage cushions



Patented massage cushions of the Philips Avent breast pump flex in and out, replicating baby's suckling action and is intended to help stimulate a fast natural let down.

## Infinitely variable control

Enjoy complete control over speed, vacuum and pattern – you can even set the pump to match your baby's own suckling action.

## Intuitive control handle

Everything that controls the pump is at your fingertips and close to the breast, which means you can make adjustments instantly and intuitively - no need to reach elsewhere.

## Patented Let-down Massage Cushion



Five, soft petal massagers gently flex in and out working in tandem with the reliable vacuum. This is intended to help stimulate a fast, natural milk flow.

## Pumps more milk faster

Double pumping has been shown to improve milk production. The twin electronic breast pump is the only breast pump designed to maintain the same level of vacuum and control for both pumps, which helps to ensure effective use.

## Quiet, effortless and portable

Quiet, effortless and portable

# Specifications

## Country of origin

England

## Development stages

Stage: 0 - 6 months

## Power

Voltage: 100 - 240 V

## What is included

Extra soft newborn flow nipple: 2 pcs  
Bottle stand/funnel cover: 4 pcs  
Breast milk container (125 ml/ 4 oz): 2 pcs  
Sealing disc for milk storage: 2 pcs  
Newborn nipple travel pack: 2 pcs  
Cool pack: 4 pcs  
Spare parts: 1 pcs

## What is included

Twin Electronic Breast pump: 1 pcs  
Microfiber travel bag: 1 pcs  
Thinsulate™ carrier: 2 pcs

