

Philips
Digital baby thermometer
set

SCH530/86



Easy and reliable temperature measurement

Digital thermometer soother for easy and accurate measurement.

Fast and reliable measurement

- Digital thermometer for easy and accurate measurement
- Quick and reliable temperature measurement within seconds

Convenience for you and your baby

- Dimpled shield prevents skin irritation
- Includes a soother in the size for a newborn baby
- Designed especially for newborn babies

Baby's safety

- Teat suitable for sterilising or boiling
- Allows boiling or sterilising without damaging electronics

PHILIPS

Specifications

Power

- Operating time on battery: 200 hour(s)
- Battery type: LR41
- Number of batteries: 2
- Removable/replaceable

Technical specifications

- Operating temperature range: 32–44 °C

Accessories

- User manual
- Battery

Technical specifications

- Accuracy: 36.5–39.0°C: $\pm 0.1^{\circ}\text{C}$
- Measurement time: 10 s

Convenience

- Automatic shut-off
- Ready indicator

Convenience

- Flexible tip
- Suitable for newborns

Highlights

Digital body thermometer

Fast, reliable and accurate digital thermometer for measuring the baby's temperature.

Measure within seconds

Never too hot, never too cold — the "Just Right" digital thermometer quickly tells you if the milk or food is safe to eat or drink. It takes the guessing out of warming bottles or jars and is great for checking if food is still warm enough after taking a short break in feeding.

Dimpled shield on soother

Dimples make sure air flows behind the teat to protect your baby's sensitive skin from irritation.

Silicon teat

The soother's specially shaped silicon teat is designed to fit your baby's mouth and will not inhibit natural development. Use it with confidence to pacify your baby.

Soother included

You can put the soother regularly in your baby's mouth to stimulate him or her to get used to it. It works like any normal soother.

Detachable electronics

By detaching the sensitive temperature measuring electronics from the soother, you can safely sterilise or boil the part your baby comes into contact with.

Sterilise teat

Thanks to the special material, you can safely sterilise the soother's teat in a microwave oven or electrical steriliser or boil it without harming the soother.



Issue date 2021-02-15

Version: 5.0.1

© 2021 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com