


PHILIPS

AVENT

Twin electronic
breast pump

SCF304/02

 Product Picture

Designed for comfort

Inspired by nature

Our unique BPA-Free pump is designed with an electronic memory that puts you in control by learning and continuing your personal pumping style.

Benefits and features:

- Ultra comfortable and effective
- Control at your fingertips
- Infinitely variable speed, vacuum and rhythm – no pre-sets
- Quiet, effortless and portable
- Pumps more milk faster

Stimulate and Express

- Gentle vacuum mimics baby's suckling for steady milk flow
- Soft massage cushions trigger natural let-down, like baby

Highlights

Soft massage cushions



Patented massage cushions of the Philips Avent breast pump flex in and out, replicating baby's suckling action and is intended to help stimulate a fast natural let down.

Specifications

Power

Voltage: 100 - 240 V

Country of origin

England

What is included

Bottle stand/funnel cover: 4 pcs

Extra soft newborn flow nipple: 2 pcs

Breast milk container (125 ml/ 4 oz): 2 pcs

Sealing disc for milk storage: 2 pcs

Twin Electronic Breast pump: 1 pcs

Cool pack: 4 pcs

Newborn nipple travel pack: 2 pcs

Spare parts: 1 pcs

Microfiber travel bag: 1 pcs

Thinsulate™ carrier: 2 pcs

Development stages

Stage: 0 - 6 months

