


PHILIPS

InfraCare infrared
lamp

650W

Half body treatment

63.5cm high



HP3641



Effective pain relief

Deep penetrating infrared warmth

Philips InfraCare gives effective pain relief for muscles and joints. Its 650 W infrared warmth is very comfortable and penetrates deep into the skin, stimulates blood circulation and gives effective relief for half-body areas of 60x40 cm.

Easy to use

- Adjustable angle with rotation

Effective pain relief

- Deep penetrating Infrared warmth

Half-body treatment

- 650 Watt infrared halogen lamp


Highlights

650 Watt infrared halogen lamp



The Philips InfraCare HP 3643 has been optimised for effective relief of pains in areas covering 60 x 40 cm, such as the entire back, both shoulders and neck or thigh. The InfraCare features innovative infrared halogen lamp technology. The special optics, filter and the powerful halogen lamp of 650 Watt have been optimised, enabling half-body treatment (60 x 40cm), and much more comfortable warmth, as the infrared warmth is more equally divided over the treatment area.

Adjustable angle with rotation



The InfraCare HP 3643 is easy to use for half-body-treatment and for different situations (e.g. lying on a bed or couch) due its maximum

flexibility. Next to the height adjustment, the lamphouse can easily be positioned to optimally treat the right body part, by moving the lamphouse into a horizontal or vertical position (left/right 90 degrees) and by tilting it up or down (45 degrees).

Deep penetrating Infrared warmth



Infrared light has been proven to effectively relieve aches and pains caused by muscular pains and stiff joints. The comfortable warmth created by infrared light penetrates deep into your skin, stimulates blood circulation and warms your muscles. As your muscles are soothed, they automatically loosen up and relax. Because it makes the tissue more pliable, warmth also reduces stiffness and makes joints more flexible.

Specifications

Easy to use

2 handgrips: for easy transport and positioning

On/off switch

Digital user interface: for feedback on time and warmth intensity

Extendable height: 45-80 cm

Easy positionable lamphouse (left/right): 90

Easy positionable lamphouse (up/down): 45

Weight and dimensions

Product weight: 5.5 kg

Product dimensions: 29.0x63.5x29.0 (WxHxD) cm

Technical specifications

Power: 650 W

Voltage: 230 V

Cord length: 300 m

Frequency: 50/60 Hz

Weight and dimensions

Qty. on Euro pallet: 18 pcs

F-box dimensions (=A-box): 33.9x80.0x36.5 (WxHxD) cm

F-box weight (=A-box): 6.5 kg

Technical specifications

Insulation: Class II (double isolation)

Lifetime of lamps: 500 hour(s)

Half-body treatment

Infrared halogen lamp: 650 W

Treatment area: 60x40 cm

Safety

IEC certified: Complies to IEC 60601 and 60335

Logistic data

Country of origin: Germany

Medical appliance

Medical Device Directive: 2007/47/EC, MDD 93/42/EEC

Logistic data

CTV code: 884364101000

Effective pain relief

Deep penetrating warmth

Treatment control

Digital timer: with auto shut-off

Dimmer with 3 warmth settings

Pause function

