

PHILIPS

Infrared lamp



HP3614



Relieves muscular pains and common colds

Deep, focused warmth

Philips infrared lamp gives pain relief for muscles and joints. Its 100 Watt infrared warmth is comfortable and penetrates deep into skin, improving blood circulation thereby giving relief for focused areas of 20x30 cm.

Philips - the light expert

- More than 15 million infrared lamps sold

Highlights

Philips - the light expert

Philips is one of the global leading companies in light, with over 100 years of light experience. Globally Philips already sold more than 15 million of Infrared lamps.

Specifications

Easy to use

Adjustability: 0-40 degrees (backwards)

Weight and dimensions

Product weight: 0.5 kg

Product dimensions: 20.0x12.0x16.0 (HxWxD) cm

Technical specifications

Power: 100 W

Voltage: 220/230 V

Cord length: 1.6 m

Frequency: 50 Hz

Weight and dimensions

A-box dimensions: 22.4x26.2x58.8 (HxWxD) cm

F-box weight: 1 kg

Qty. on Euro pallet: 288 pcs

No. of F-boxes in A-box: 6 pcs

A-box weight: 3.1 kg

F-box dimensions: 21.1x12.6x19.2 (HxWxD) cm

Technical specifications

Insulation: Class II (double isolation)

Type of lamps: R95, 100W

Lifetime of lamp: 750 sessions of 10 minutes

Safety

IEC certified: In line with IEC 60335

Logistic data

Country of origin: Germany

Logistic data

CTV code: 884361401000

Easy storage

Cord storage

