

Philips AVENT  
Digital baby thermometer  
set



SCH540/00

## Gentle and accurate

Combining a digital and a soother baby thermometer

Helping you to get a measurement of your baby's temperature, without upsetting it. The soother thermometer SCH540/00 uses the same shield and teat as the AVENT Freeflow soother for easy transition from soother to temperature measurement.

### **Convenient for you and your baby**

- Includes AVENT Freeflow Pacifier

### **Fast and reliable digital temperature readings**

- Clear measurement displays and battery status indicators
- Battery low indicator

### **Comfort for your child**

- Soft, flexible tip digital thermometer
- Choice to adapt temperature measurement method to suit baby
- Specially designed pacifier thermometer

**PHILIPS**  
**AVENT**

# Specifications

## Country of origin

- Made in China: Thermometers, Yes

## Development stages

- Stage: 0 - 6 months, 6 - 12 months

## Power

- Operating time on battery: >100 hour(s)
- Battery type: LR41
- Number of batteries: 1 per thermometer

## Technical specifications

- Operating temperature range: 32-42.9 °C

## Accessories

- User manual



## Technical specifications

- Accuracy: 35-42 +/- 0.1 °C
- Measurement time: 10secs ±3secs in water (bath), 25secs ±10secs oral & rectal, 60secs ±10secs underarm

## Convenience

- Automatic shut-off

## Country of origin

- Made in England: Soothers, Yes

## Convenience

- Flexible tip
- Suitable for newborns
- Clear measurement display
- Battery status indicator

# Highlights


## Clear measurement displays

Clear measurement displays and battery status indicators

## Avent Freeflow Pacifier

Lets your baby get used to the shape.

## Battery low indicator



Battery low indicator

## Choice of measurement method

The Philips Avent thermometer accurately and conveniently measures your baby's body temperature orally, axially and rectally.

## Soft, flexible tip digital thermometer

Soft, flexible tip digital thermometer

## Specially designed pacifier thermometer

Enables easy transition from pacifier to temperature measurement.



Issue date 2021-02-15

Version: 9.1.1

© 2021 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)