



Philips  
InfraPhil infrared lamp

150W



HP3616

## Relieves muscular pains

Deep, focused infrared warmth

Deep, focused warmth from infrared lamp

### Easy to use

- On/off switch
- Adjustable angle

### Relieves muscle pain

- Deep penetrating Infrared warmth

### Easy to store

- Cord storage

# PHILIPS

# Highlights

## Deep penetrating Infrared warmth



Infrared light has been proven to effectively relieve aches and pains caused by muscular pains and stiff joints. The comfortable warmth created by infrared light penetrates deep into your skin, stimulates blood circulation and warms your muscles. As your muscles are soothed, they automatically loosen up and relax. Because it makes the tissue more pliable, warmth also reduces stiffness and makes joints more flexible.

## Adjustable angle



The infrared lamp can easily be positioned up to 40 degrees backwards.

## On/off switch



Easy on/off switch at the back of the infrared lamp, so the appliance does not have to be unplugged after each treatment.

## Cord storage



The infrared lamp has a cord storage holder at the back of the appliance. Simply wrap the cord around it and it can easily be stored again.

# Specifications

## Easy to use

- Adjustability: 0-40 degrees (backwards)
- On/off switch

## Weight and dimensions

- Product dimensions: 22.0x13.0x18.5 (HxWxD) cm
- Product weight: 1 kg

## Technical specifications

- Voltage: 220/230 or 240/250V V
- Power: 150 W
- Cord length: 180 m
- Frequency: 50 Hz

## Weight and dimensions

- A-box dimensions: 24.0x28.8x62.4 (HxWxD) cm
- F-box weight: 1 kg
- A-box weight: 6.4 kg
- Qnt. on Euro pallet: 192 pcs
- No. of F-boxes in A-box: 6 pcs
- F-box dimensions: 22.5x14.0x20.4 (HxWxD) cm

## Technical specifications

- Insulation: Class II (double isolation)
- Type of lamps: PAR 38E, 150 W + prismatic rings for more focus
- Lifetime of lamp: 750 sessions of 10 minutes

## Safety

- IEC certified: Complies to IEC 60335

## Logistic data

- Country of origin: Germany

## Medical appliance

- Medical Device Directive: 2007/47/EC, MDD 93/42/EEC

## Logistic data

- CTV code: 884361601000

## Easy storage

- Cord storage

