

# PHILIPS

## AVENT

Twin electronic  
breast pump

SCF304/60

 Product Picture

## Designed for comfort

Inspired by nature

Our unique BPA-Free pump is designed with an electronic memory that puts you in control by learning and continuing your personal pumping style.

### **Benefits and features:**

- Learns and continues your pumping rhythm
- Ultra comfortable and effective
- Infinitely variable speed, vacuum and rhythm – no pre-sets
- Quiet, effortless and portable
- Pumps more milk faster

### **Stimulate and Express**

- Soft massage cushions trigger natural let-down, like baby
- Gentle vacuum mimics baby's suckling for steady milk flow

# Highlights

**Learns and continues your pumping rhythm**

Electronic memory learns and at the touch of a button, continues your personal pumping rhythm.

**Gentle-draw vacuum**

The gentle-draw vacuum mimics baby's suckling to create a steady milk flow that requires less pumping.

**Infinitely variable control**

Enjoy complete control over speed, vacuum and pattern – you can even set the pump to match your baby's own suckling action.

**Patented Let-down Massage Cushion**

Five, soft petal massagers gently flex in and out working in tandem with the reliable vacuum. This is intended to help stimulate a fast, natural milk flow.

**Pumps more milk faster**

Double pumping has been shown to improve milk production. The twin electronic breast pump is the only breast pump designed to maintain the same level of vacuum and control for both pumps, which helps to ensure effective use.

**Quiet, effortless and portable**

Quiet, effortless and portable

**Soft massage cushions**

Patented massage cushions of the Philips Avent breast pump flex in and out, replicating baby's suckling action and is intended to help stimulate a fast natural let down.

# Specifications

**Country of origin**

England

**Development stages**

Stage: 0 – 6 months

**Power**

Voltage: 100 – 240 V

**What is included**

Extra soft newborn flow nipple: 2 pcs  
Bottle stand/funnel cover: 4 pcs  
Sealing disc for milk storage: 2 pcs  
Cool pack: 2 pcs  
Breast milk container (260 ml/ 9 oz): 2 pcs

**What is included**

Twin Electronic Breast pump: 1 pcs  
Manual pump parts for when power is not available: 1 pcs  
Microfiber travel bag: 1 pcs  
Thinsulate™ carrier: 2 pcs

