

PHILIPS

 Product Picture

AVENT

Twin electronic
breast pump

SCF294/02

Designed for Comfort

Inspired by nature

Our unique pump is designed with an electronic memory that puts you in control by learning and continuing your personal pumping style.

Benefits and features:

- Ultra comfortable and effective
- Control at your fingertips
- Infinitely variable speed, vacuum and rhythm – no pre-sets
- Quiet, effortless and portable
- Pumps more milk faster

Stimulate and Express

- Soft massage cushions trigger natural let-down, like baby
- Gentle vacuum mimics baby's suckling for steady milk flow

Highlights

Soft massage cushions



Patented massage cushions of the Philips Avent breast pump flex in and out, replicating baby's suckling action and is intended to help stimulate a fast natural let down.

Gentle-draw vacuum

The gentle-draw vacuum mimics baby's suckling to create a steady milk flow that requires less pumping.

Patented Let-down Massage Cushion



Five, soft petal massagers gently flex in and out working in tandem with the reliable vacuum. This is intended to help stimulate a fast, natural milk flow.

Pumps more milk faster

Double pumping has been shown to improve milk production. The twin electronic breast pump is the only breast pump designed to maintain the same level of vacuum and control for both pumps, which helps to ensure effective use.

Infinitely variable control

Enjoy complete control over speed, vacuum and pattern – you can even set the pump to match your baby's own suckling action.

Intuitive control handle

Everything that controls the pump is at your fingertips and close to the breast, which means you can make adjustments instantly and intuitively- no need to reach elsewhere.

Quiet, effortless and portable

Quiet, effortless and portable

Specifications

Country of origin

England

Development stages

Stage: 0 - 6 months

Power

Voltage: 100 - 240 V

What is included

- Extra soft newborn flow nipple: 2 pcs
- Breast milk container (125 ml/ 4 oz): 2 pcs
- Bottle stand/funnel cover: 4 pcs
- Sealing disc for milk storage: 2 pcs
- Newborn nipple travel pack: 2 pcs
- Cool pack: 4 pcs
- Spare parts: 1 pcs
- Breast milk container (260 ml/ 9 oz): 2 pcs

What is included

- Twin Electronic Breast pump: 1 pcs
- Manual pump parts for when power is not available: 1
- Microfiber travel bag: 1 pcs
- Thinsulate™ carrier: 2 pcs

