

Philips Avent Express
Electric Steam Steriliser

220-240 V

SCF274/23



Sterilises 6 bottles in 8 minutes

Fast and easy to use

Fast and easy to use

Benefits and features:

- Fast and easy to use
- Sterile contents for up to 6 hours
- Effective sterilisation
- Large capacity

PHILIPS
AVENT

Electric Steam Steriliser
220-240 V

Specifications

Country of origin

- England

Development stages

- Stage: 0–6 months

Power

- Voltage: 220 - 240 V

What is included

- Extra soft Newborn Flow teat: 3 pcs
- Extra soft Slow Flow teat: 1 pcs

Weight and dimensions

- Dimensions: 313 (H), 235 (W), 235 (L) mm
- Weight: 1.575 kg

What is included

- Tongs: 1 pcs
- Bottle and teat brush: 1 pcs
- Feeding Bottle (260 ml/9 oz): 1 pcs
- Feeding Bottle (125 ml/4 oz): 1 pcs
- Newborn Soother: 2 pcs

What is included

- Measuring jug: 1 pcs
- Express Electric Steam Steriliser: 1 pcs

SCF274/23

Highlights

Fast and easy to use



Just add water, load and switch on. Contents are sterile and ready to use in approximately 8 minutes.

Sterile contents for up to 6 hours



Contents remain sterile for up to 6 hours if unopened.

Large capacity



Holds up to six 260 ml/9 oz Avent Bottles or two Philips Avent Breast Pumps and accessories.

Effective sterilisation



Steam sterilisation is proven to be the most effective way to protect your baby from harmful germs. Based on the hospital principle, the intensive heat of the steam eliminates harmful bacteria.



Issue date 2021-02-12

Version: 6.3.1

© 2021 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com