

PHILIPS

IPL hair removal
system

Lumea Essential

For use on body only

15 minutes to treat lower legs

Lifetime >100.000 light
pulses

Extra long cord



SC1991/00



The ideal way to prevent hair regrowth for smooth skin at home

Philips Lumea IPL works wonders to prevent the reappearance of hair on body. Gentle pulses of light, applied regularly, keep skin silky-smooth every day.

Long lasting smoothness

- Effortlessly effective
- Enjoy smooth skin every day

Hair removal with IPL technology

- Adapted for safe use in the comfort of your own home

Gentle and effective, even on sensitive skin

- For use on the body
- Adjustable to your skin type
- Five adjustable light energy settings
- Bigger body attachment for fast treatment

Compact corded design

- No replacement parts, no hidden costs

Highlights

Enjoy smooth skin every day



Philips Lumea applies gentle pulses of light to the hair root. This results in the hair shedding naturally and hair regrowth is inhibited. Repeating this treatment regularly will keep your skin smooth everyday.

Effortlessly effective



Our clinical studies have shown significant hair reduction in as little as four bi-weekly treatments resulting with smooth hair free skin. To maintain these results, simply repeat the treatment when needed. Time between treatments may vary based on your individual hair regrowth.

IPL technology



Philips Lumea uses an innovative light-based technology called IPL (Intense Pulsed Light), derived from the technology used in

professional beauty salons. Philips has adapted this technology for safe and effective use in the comfort of your own home. Philips worked closely with leading dermatologists to develop our breakthrough hair removal system. For more than 10 years we carried out extensive consumer research with more than 2000 volunteers.

5 light energy settings



Philips Lumea has five adjustable light energy settings to ensure gentle but effective treatment. Used correctly, Philips Lumea's IPL treatment is safe and gentle to use even on sensitive skin and sensitive body areas.

Application areas



Use with confidence to remove hair on body areas below the neckline: legs, armpits, bikini area, arms and belly.

Bigger body attachment, 4cm2




Bigger body attachment enables fast treatment on large areas like legs.

No replacement parts needed



Philips Lumea is a complete solution that comes ready to work straight out of the box. More importantly, it does not require any replacement lamps or gels.

Suitable hair and skin types



Philips Lumea works effectively on (naturally) dark blonde, brown and black hairs. As with other light-based treatments, Philips Lumea is NOT effective on red, light-blonde or white/grey hair. Philips Lumea is also NOT suitable for dark skin.

Specifications

Application areas

Body areas: Arms, Legs, Belly, Bikini, Armpits

Safety and adjustable settings

Integrated UV filter: Protects skin from UV light

5 light energy settings: Adjustable to your skin type

Integrated safety system: Prevent unintentional flashing

Application mode

Step and Flash: For treatment on curvy areas

Corded/cordless use: Corded

Power

Voltage: 100-240 V

Technical specifications

High performance lamp: Delivers <gt/> 100,000 flashes

Attachments

Body attachment (4 cm²): For use below the neckline

Service

Warranty: 2 years global warranty

Items included

Instructions for use: Quick Start Guide, User manual

Storage: Pouch

Adapter: 24V / 1500mA

Application time

Lower leg: 8 min

Armpit: 1 min

Bikini line: 1 min

