

PHILIPS

AVENT

3-in-1 electric steam steriliser

Kills 99.9% of harmful germs

Sterilises in 6 minutes

Fits 6 Philips Avent bottles

Adjustable 3-in-1 design

SCF285/01



Convenient and effective sterilisation

Flexible, easy loading

With its adjustable size, the Philips Avent 3-in-1 electric steam steriliser takes up the least amount of kitchen space, yet perfectly fits the items you want to sterilise, whether it's a few small items or a full load.

Effective sterilising using natural steam

- Sterile for up to 24 hours if the lid is unopened

Fits various bottle types and other products

- Sterilises various bottles, breast pumps and accessories
- Spacious design

Space-saving 3-in-1 design

- 3-in-1 modular steriliser design

Highlights

Flexible 3-in-1 steriliser



The steriliser's unique modular design makes it quick and easy to fit and organise your bottles and accessories. It also takes up minimal kitchen space.

Sterile for up to 24 hours



The steriliser will keep its contents—baby bottles, breast pumps, etc.—sterile for up to 24 hours if the lid is unopened.

Spacious design



Holds up to six 330 ml / 11 oz Philips Avent Classic and Natural bottles to be sterilised at once.

Works across bottle types

The steriliser is suitable for use with both standard-neck and wide-neck baby bottles. It also fits other baby products like breast pumps and accessories.

Power

Voltage: 220-240 V

Technical specifications

Voltage: 50-60Hz

Power consumption: 650 W

Safety Classification: Class 1

Weight and dimensions

Dimensions: 290 x 160 x 150 (w x d x h) mm

Weight: 1.5 kg

Country of origin

Made in: Turkey

Material

Polypropylene

What is included

Newborn soother: 1 pcs

Tongs: 1 pcs

Electric steam steriliser: 1 pcs

Classic Bottle 4 oz: 1 pcs

Compatibility

Philips Avent range-compatible

Development stages

Stage: 0-6 months

