



Philips Avent
Large massage cushion

Fits all comfort breast pump



SCF167/01

Right size, maximum comfort

Breast pump with wider funnel

This soft massage cushion has a wider funnel, specially designed for mums with larger nipples. It can be used with Philips Avent Comfort breast pump range. See below images to help you decide if you need the larger breast pump cushion.

Comfortable pumping

- Wider funnel for comfortable fit
- Fits all Natural range breast pumps

Check your fit for optimal comfort

- Right fit: The nipple sits comfortably in the cushion funnel
- Choose large size: Nipple touches the walls of the funnel

Large massage cushion
Fits all comfort breast pump

SCF167/01

Specifications

Country of origin

- England

Development stages

- Stage: Pregnancy, 0 - 6 months

Material

- BPA free*

Replaceable part

- Fits product types :: SCF330, SCF332, SCF334, SCF330, SCF332

What is included

- Large size cushion: 1 pcs

Functions

- Soft massage cushion: Gentle stimulation

Cushion

- Material: Silicone, BPA free*

Compatibility

- Compatible with Natural range

Ease of use

- Cushion: Easy to clean, Dishwasher proof


Highlights

Large breast pump cushion



You can find out whether you need a large breast pump cushion by checking whether your nipple fits comfortably in the cushion funnel. If your nipple touches or rubs against the walls of the cushion funnel, then you should choose the large breast pump cushion to ensure a more comfortable expressing experience.

Correct fit



If the nipple sits comfortably in the cushion funnel, you are using the correct size massage cushion.

Too small fit



If your nipple touches or rubs against the wall of the cushion funnel, then you should choose the large breast pump cushion with wider funnel to ensure a more comfortable pumping experience.

Fits all Natural breast pumps

The large size soft massage cushion is compatible with all Philips Avent Comfort breast pumps in the Natural range.



Issue date 2021-03-09

Version: 6.0.1

EAN: 08 71010 35779 73

© 2021 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com