

PHILIPS

AVENT

Manual Breast Pump and Store Set

Express, store and feed.

Pump with soft massage cushion

Includes 2 bottles and teats

10 x storage cups and lids

SCD221/00



More comfort, more milk* Pump, store and feed

Manual breast pump set with storage cups

The Comfort Manual Breast Pump & Store set is ideal for mums who plan to express their breast milk. The convenient set includes all you need to pump and store your breast milk and feed your baby, and includes 2 bottles and 10 storage cups.

Easy to combine breast and bottle feeding

- Includes our Natural bottle and teat for natural latch-on

Easy to set up, express and clean

- Few separate parts and intuitive design
- Compact lightweight design

Clinically proven less leaning forward*

- Relaxing position thanks to its special design

Store - Storage cups

- For safe storage and transport
- Pump in storage containers

Soft massage cushion with massaging petals

- Soft massage cushion gently stimulates milk flow

Highlights

Relaxing position



The breast pump has a unique design, so your milk flows directly from your breast into the bottle, even when you are sitting up straight. This means you can sit more comfortably when pumping; no need to lean forward to make sure all your milk ends up in the bottle. Sitting comfortably and being relaxed when pumping helps your milk to flow more easily.

Soft massage cushion



Our massage cushion has a soft, velvety texture that gives a warm feel to the skin for comfortable, gentle stimulation of your milk flow. The iconic petal cushion is designed to mimic baby suckling, intended to gently help stimulate let-down.

Easy on-the-go



The breast pump is small and lightweight, which means it is easy to store and transport, making pumping on the go more discreet.

Easy to set up, express and clean



Few separate parts and intuitive design. Easily express milk for your baby, anytime.

Includes our Natural bottle



The wide, breast-shaped teat promotes natural latch-on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

Pump in storage containers



With the milk storage cup adapter, you can easily pump in the milk storage container. No leakage, no spillage of your precious breast milk.

Leak-proof lid



Philips Avent storage cups have a lid that creates a secure seal for safe storage and transport.

Specifications

Design

Bottle design: Ergonomic shape, Wide neck
Breast pump design: Compact design

Material

Bottle: BPA-free*, Polypropylene
Teat: BPA-free*, Silicone
Breast pump: BPA-free* (food contact parts only)
Breast pads: Natural breathable material, Silky soft top-sheet
Storage cups: BPA free*, Polypropylene

What is included

Breast pump body with handle: 1 pcs
Milk storage cup adapter: 2 pcs
Sealing disc for milk storage: 1 pcs
Natural bottle 125 ml/4 oz: 1 pcs
Natural bottle 9 oz: 1 pcs
Milk storage cup (180 ml/6 oz): 10 pcs
Breast pad sample packs: 1 pack (2 pcs)
Standard size cushion (19.5 mm): 1 pcs (A larger-size cushion is available and sold separately)

Ease of use

Breast pump use: Easy cleaning, Fully compatible range, Intuitive assembly

Functions

No leaning forward: Sit in a comfortable position
Soft massage cushion: Gentle stimulation

Development stages

Stage: 0–6 months



* Clinically proven comfort: Tests carried out on 110 mothers showed a significant preference for Philips Avent vs. a leading competitor
* More comfort: 73% of 73 UK breastfeeding mothers agreed that this breast pump was more comfortable than their current pump (main market brands).
* Independent research has shown that there may be a link between stress levels and milk production. See www.philips.com/Avent
* BPA-free breast pump: Only associated with the bottle and other parts that come into contact with breast milk. Following EU regulation, 10/2011