

Philips Philips Avent
Electric breast pump
diaphragm

Breast pumps

CP9914/01



Creates the vacuum needed for expression

Check for compatibility below

This silicone diaphragm creates the vacuum needed for expressing breast milk. For use with electric breast pumps.

Ease of using Philips spare parts

- Easily renew your product with original Philips parts

PHILIPS

Specifications

Replaceable part

- Fits product types: SCD292/31 SCD395/21, SCF332/01 SCF332/01 SCF332/11 SCF332/31 SCF332/33 SCF332/61 SCF332/62, SCF334/02 SCF334/03 SCF334/04 SCF334/12 SCF334/14,

SCF334/31 SCF334/61 SCF391/11 SCF391/71 SCF392/11, SCF393/11 SCF393/61 SCF393/82 SCF394/11 SCF394/52 SCF394/61, SCF395/11 SCF395/31 SCF395/91, SCF396/11 SCF396/31, SCF397/11 SCF398/11 SCF398/31

Highlights

Easy replacement of product parts

From time to time your product needs a facelift, and with Philips consumer replacement parts to renew your product, it's never been so easy! All this with guaranteed Philips quality.



Issue date 2024-04-13

Version: 5.5.1

EAN: 87 10103 62314 4

© 2024 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com