



Philips Avent  
Comfort Single electric  
breast pump

### Natural

Includes 4oz bottle



SCF332

## More comfort, more milk naturally

### Breast pump with soft massage cushion

When you are comfortable and relaxed, your milk flows more easily. That is why we created our most comfortable breast pump yet: Sit comfortably with no need to lean forward and let our soft massage cushion gently stimulate your milk flow.

#### Sit more comfortably with no need to lean forward

- More comfortable pumping position due to unique design

#### Easy to combine breast and bottle feeding

- Includes our Natural bottle and nipple for natural latch on

#### Gently stimulates let down and milk flow

- Soft massage cushion with warm feel

#### Other benefits

- Intuitive assembly. Easy visual matching of parts
- Compact lightweight design
- Compatible with other Philips Avent feeding products
- Easy cleaning due to the small number of separate parts
- Easy operation at the touch of a button

#### Choose your most comfortable setting

- Features a gentle stimulation mode and 3 pumping settings

**PHILIPS**  
**AVENT**

Comfort Single electric breast pump  
Natural Includes 4oz bottle

SCF332/60

# Highlights

## More comfortable position



The breast pump has a unique design, so your milk flows directly from your breast into the bottle, even when you are sitting up straight. This means you can sit more comfortably when pumping; no need for you to lean forward to make sure all your milk ends up in the bottle. Sitting comfortably and being relaxed when pumping, naturally helps your milk to flow more easily.

## Simple settings to choose from



When switched on, the pump automatically starts in gentle stimulation mode to get your milk flowing. Then choose from 3 pumping settings to make milk flow most comfortable for you.

## Soft massage cushion



Our massage cushion has a new soft velvety texture that gives a warm feel to the skin for comfortable, gentle stimulation of your milk flow. The cushion is designed to gently mimic your baby's suckling to help stimulate let down.

## Includes our Natural nipple



The wide breast shaped nipple promotes natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

## Compact design



The breast pump has a compact design, which makes it easy to hold and position on your breast. The small, lightweight base unit can easily be placed within comfortable reach for full control when pumping. For extra

convenience in transport and storage, the tube simply wraps around the base unit. Easy to use on the go with batteries.

## Easy cleaning

Cleaning is easy, thanks to the small number of separate parts. Your milk will never come into contact with the tubing and base unit. All parts are dishwasher proof, except the electrical parts.

## Easy operation

A single electric breast pump is ideal for moms who pump milk frequently and value being able to pump effortlessly, at the touch of a button.

## Intuitive assembly



Easy visual matching of parts for intuitive assembly.

## Fully compatible range



This breast pump can be used in combination with other feeding products in the Philips Avent range, including our Classic bottles and milk storage containers. Philips Avent also provides a range of breast care accessories to help you breastfeed for longer and enhance your comfort.

# Specifications

## Country of origin

- England

## Development stages

- Stage: 0 - 6 months

## Material

- Breast pump: Polypropylene, BPA free\*
- Bottle: Polypropylene, BPA free\*
- Nipple: Silicone, BPA free\*

## What is included

- Extra soft newborn flow nipple: 1 pcs
- Sealing disc for milk storage: 1 pcs

## Design

- Breast pump design: Compact design
- Bottle design: Ergonomic shape, Wide neck

## What is included

- Breast pad sample packs: 2 (1 Day & 1 Night

- pack) pcs
- Natural bottle 4oz: 1 pcs
- Breast pump body: 1 pcs
- Standard size cushion: 1 pcs
- Base unit incl. tubing: 1 pcs
- Spare diaphragm: 1 pcs
- Travel cover: 1 pcs

## Ease of use

- Breast pump use: Intuitive assembly, Easy cleaning, Fully compatible range

## Breast pump

- Material: BPA free\*

## Functions

- No leaning forward: Sit in a comfortable position
- Soft massage cushion: Gentle stimulation
- Settings: 3 Expression settings, 1 Stimulation mode



Issue date 2021-03-10

Version: 4.1.1

EAN: 08 71010 36169 48

© 2021 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)

\* 0% BPA, following EU regulation 10/2011