

# PHILIPS

## AVENT

### Manual Breast Pump

Ultra Comfort

Includes 125 ml (4 oz) Natural bottle



SCF330/70



## More comfort, more milk\*

### Clinically proven comfort\*

When you are comfortable and relaxed, your milk flows more easily. That is why we created our most comfortable breast pump yet: sit comfortably with no need to lean forward and let our soft massage cushion gently stimulate your milk flow.

#### **Easy to combine breast and bottle feeding**

- Includes our Natural bottle and teat for natural latch-on

#### **Easy to set up, express and clean**

- Few separate parts and intuitive design
- Compact lightweight design

#### **Clinically proven less leaning forward\***

- Relaxing position thanks to its special design

#### **Soft massage cushion with massaging petals**

- Soft massage cushion gently stimulates milk flow

# Highlights

## Relaxing position



The breast pump has a unique design, so your milk flows directly from your breast into the bottle, even when you are sitting up straight. This means you can sit more comfortably when pumping; no need to lean forward to make sure all your milk ends up in the bottle. Sitting comfortably and being relaxed when pumping helps your milk to flow more easily.

## Soft massage cushion



Our massage cushion has a soft, velvety texture that gives a warm feel to the skin for comfortable, gentle stimulation of your milk flow. The iconic petal cushion is designed to mimic baby suckling, intended to gently help stimulate let-down.

## Easy on-the-go



The breast pump is small and lightweight, which means it is easy to store and transport, making pumping on the go more discreet.

## Easy to set up, express and clean



Few separate parts and intuitive design. Easily express milk for your baby, anytime.

## Includes our Natural bottle



The wide, breast-shaped teat promotes natural latch-on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

# Specifications

## Design

**Bottle design:** Ergonomic shape, Wide neck  
**Breast pump design:** Compact design

## Material

**Bottle:** Polypropylene  
**Teat:** Silicone

## What is included

**Extra-soft newborn flow teat:** 1 pcs  
**Breast pump body with handle:** 1 pcs  
**Sealing disc for milk storage:** 1 pcs  
**Natural bottle 125 ml/4 oz:** 1 pcs  
**Travel cover:** 1 pcs  
**Standard massage cushion:** 1 pc (19.5 mm [0.77"], a larger size cushion is available separately)  
**Disposable breast pads:** 12 pcs

## Ease of use

**Breast pump use:** Easy cleaning, Fully compatible range, Intuitive assembly

## Functions

**No leaning forward:** Sit in a comfortable position  
**Soft massage cushion:** Gentle stimulation

## Development stages

**Stage:** 0–6 months

