

PHILIPS

AVENT

Silicone diaphragm
for breast pumps

CP9914



To bring vacuum in the breast pump.

Find matching products on the specifications tab

This silicone diaphragm creates vacuum in your electric breast pump.

Ease of using Philips spare parts

- Easy renew your product with original Philips parts

Highlights

Easy renew of product parts

From time to time your product needs a facelift but it was never so easy like with Philips consumer replacement parts to renew your product. This all with guaranteed Philips quality.

Specifications

Replaceable part

Fits product types :: SCF394, SCF398, SCF396, SCF397, SCF332, SCD292/31, SCF393, SCF395, SCF391, SCF332, SCF334

