



Philips Avent
DECT baby monitor



SCD530/00

Healthy climate, total comfort

Zero interference guaranteed with DECT technology

Enjoy complete peace of mind with the new Philips Avent range of DECT Baby Monitors you can relax, safe in the knowledge your baby's happy even when you are out the room.

Zero interference

- DECT technology guarantees zero interference

Crystal clear sound

- LEDs light to signal noise levels
- Choose the best sound level for your surroundings
- Reassurance with crystal clear sound

330 m range

- Up to 24 hour operation time before recharge
- Always know that you are connected
- Extra long range for total mobility

Complete comfort

- Stay connected with talk back
- Remote activated night light and lullabies
- Monitor the humidity & temperature in the baby's room

Highlights

Guaranteed zero interference



DECT technology guarantees zero interference from any other transmitting product like other baby monitors, cordless phones and cell phones. Data encryption provides a secure and private connection, so you can be sure that you are the only one that can hear your baby

Reassuring crystal clear sound



Hear every tiny giggle, gurgle and hic-up with perfect clarity. DECT technology provides high quality, crystal clear sound so you can hear your baby anytime

Superior operating time



The small, rechargeable parent unit will give you the freedom of cordless roaming for up to 24 hours before re-charging is required

Humidity & temperature sensor




Combined humidity and temperature sensor allows you to monitor the climate in baby's room. Your baby can't regulate its body temperature as well as you can and a slight change in temperature, and too much or too little humidity can make your baby restless. A specially combined humidity and temperature sensor will immediately alert you, through the digital display on the parent unit, if the climate in your baby's room changes. Monitor the complete climate in baby's room

LEDs light to signal noise



Lights are activated to indicate the level of sound in your baby's room. Even with the Parent Unit is muted, a vibrating sound alert will notify you when a set sound level is reached

Extended operating range



With a range of 330m, a belt clip and neck cord so you can be fully mobile in and around the home. A versatile travel bag will protect your baby monitor when you are away from home

Adjustable sound sensitivity



Enables you to adjust the sound and volume levels, so you can hear every giggle, gurgle and hic-up

Specifications

Development stages

- Stage: Pregnancy, 0 - 6 months

Power

- Operating time on battery: 24 hour(s)
- Power supply: 220 - 240 V
- Charging time: 3 hour(s)

Tuner/Reception/Transmission

- Automatic channel selection
- Number of channels: 120

Convenience

- Auto out-of-range warning
- Power on indication
- Battery low indication
- Battery charging indication
- Volume control
- Sound-level lights

Technical specifications

- Operating temperature range: 10 - 40 °C

Power

- Power supply: 120 V (US)

Tuner/Reception/Transmission

- Frequency band: DECT

Accessories

- Batteries
- Charging base for parent unit
- Quick start guide

Convenience

- Sensitivity control

Accessories

- Belt clip, neck cord and travel bag
- AC/DC adaptor
- DFU/ user manual

Logistic data

- F-box dimensions (W x H x D): 305 x 135 x 187 mm

Convenience

- Vibration alert

Technical specifications

- Storage temperature range: 10 - 40 °C



Issue date 2021-03-10

Version: 12.4.1

EAN: 00 07502 00054 72

© 2021 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* This product is available in the US and Canada with only fahrenheit settings.