

PHILIPS



Under-the-nose
mask

DreamWear Full Face

Top-of-head tube design

Includes 4 point headgear design

Comfortably sits under the nose

Available in multiple sizes



HH1052/02

Like wearing nothing at all*

With the freedom and comfort of DreamWear Full Face, you can sleep how you want.

Magnetic clips

- The headgear clips and mask frame contain magnets

Open field of vision

- Open design offers a wide field of vision

Innovative design

- Freedom of movement
- A modular design for more options

Comfort and ease of use

- Easy to use
- Soft, flexible frame
- Under the nose full face cushion

DreamWear Full Face

Under-the-nose mask

HH1052/02

Highlights

Freedom of movement




By moving the tube location to the top of the head, DreamWear Full Face offers more freedom of movement during the night. Users stated they didn't need to choose a sleeping position when wearing DreamWear Full Face.*

No more nose bridge red marks



The under the nose innovative design of the cushion means no more red marks, discomfort, and irritation on the nose bridge.**

Unobstructed view



Open design makes it easy to wear glasses, read, watch TV, and use a computer or tablet before falling asleep.**

Like wearing nothing at all



The soft, flexible frame is comfortable to wear at night. Users have stated it feels like wearing nothing at all.*

Freedom of choice



The modular design allows users to choose the cushion that works best for them, all on one frame.

Easy to use



Only a few parts for easy assembly, disassembly and cleaning.

Magnetic clips



Some medical implants may be contraindicated for use with magnetic devices and this mask may not be suitable for use. Contact your healthcare professional before you use this mask.

Specifications

Warning: Magnetic interference

Beware with pacemakers: Read IFU's before use.

Intent of use

To be used on patients $<gt;/>30\text{kg}$: with CPAP prescription

Cleaning: Headgear, tubing

Hand wash with mild detergent: Rinse thoroughly, hang to dry¹

Cleaning: Mask

Hand wash with mild detergent: Rinse thoroughly, air dry¹

Replacement

Inspect daily for wear: Replace at first sign of wear

Product details

Materials: Nylon, Polyurethane foam, Polycarbonate, Polyester, Silicone, Spandex
Warranty: 90 days

Operating pressure: 4 to 30 cmH₂O

Sound level: 17 dBA

Number of parts: 7



** Data analysis of Aug 2017 Patient Preference trial where n = 85 and prescribed masks include Resmed Airfit F10 and F20, Respiration Amara View and F&P Simplus; based on 10 days of use.

*** Design claim - the design of the mask is such that it does not contact the nasal bridge.

¹Please refer to Instructions For Use for more detailed cleaning instructions.

* The scope of supply includes a standard mask frame (size M) and four cushions for all sizes (S, M, MW, L). The sleep therapy device is not included.

* The costs of DreamWear are private and are usually not reimbursed by the statutory health insurance.