

PHILIPS

AVENT

Double Electric Breastfeeding set

Express, store and feed.

Includes 2 cushion sizes

Travel bag, handle, 4 bottles

25 Milk bags and 10 cups

SCF334/26



More milk in less time* Breastfeeding made easy

This set contains all you need to pump, store and feed your breast milk. Reduce expression time with double pumping, and sit comfortably with no need to lean forward and let our soft massage cushion gently stimulate your milk flow.

Easy to set up, express and clean

- Includes our Natural bottle and nipple for natural latch on
- Few separate parts and intuitive design

Effective expression with 3 customizable settings

- Choose the most effective and comfortable setting for you

Choose from a massaging and a power cushion

- Soft massage cushion with massaging petals
- More comfortable pumping position due to unique design
- Large power cushion for extra suction

More milk in less time*

- Reduce expression time with double pumping

Highlights

More milk in less time



Reduce expression time with double pumping compared to single pumping.

Simple settings to choose from

When switched on, the pump automatically starts in gentle stimulation mode to get your milk flowing. Then choose from 3 pumping settings to make milk flow most comfortable for you.

Soft massage cushion



Our massage cushion has a soft velvety texture that gives a warm feel to the skin for comfortable, gentle stimulation of your milk flow. The iconic petal cushion is designed to mimic baby suckling, intended to gently help stimulate let down.

Large power cushion

Choose the power cushion for the moments you need a stronger suction.

Relaxing position



The breast pump has a unique design, so your milk flows directly from your breast into the bottle or container, even when you are sitting up straight. This means you can sit more comfortably when pumping: no need for you to lean forward to make sure all your milk ends up in the bottle. Sitting comfortably and being relaxed when pumping, helps your milk to flow more easily.

Easy set up, express and clean



Few separate parts and intuitive design. Easily express milk for your baby, anytime. Closed system, no milk spillage in the tubes. Designed for hygienic expressing.

Includes our Natural bottle



The wide breast shaped nipple promotes natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

Specifications

Design

Breast pump design: Compact design
Bottle design: Ergonomic shape, Wide neck

Country of origin

England

Material

Bottle: BPA free*, Polypropylene
Nipple: BPA free*, Silicone
Breast pump: BPA free* (food contact parts only)

What is included

Extra soft newborn flow nipple: 2 pcs
Sealing disc for milk storage: 2 pcs
Natural bottle 4oz: 4 pcs
Breast pump body: 2 pcs
Travel bag: 1 pcs
Base unit incl. tubing: 1 pcs
Travel cover: 2 pcs
Breast pad sample packs: 1 pack (2pcs)
Manual handle: 1 pcs
Milk storage bags: 25 pcs
Milk storage cups: 10 pcs
Large massage cushion (25mm): 2 pcs (A smaller size cushion is available separately)
Storage cups adaptors: 2 pcs
Power cushion (27mm): 2 pcs

Ease of use

Breast pump use: Easy cleaning, Fully compatible range, Intuitive assembly

Functions

No leaning forward: Sit in a comfortable position
Soft massage cushion: Gentle stimulation
Settings: 3 Expression settings, 1 Stimulation mode
More milk in less time: The perfect time saver

Development stages

Stage: 0 - 6 months



* More milk in less time: Reduce expression time with double pumping compared to single pumping.
 * 60% of US breastfeeding mothers preferred the Philips Avent Comfort Double electric breast pump for comfort when expressing vs 29% of mothers who preferred the Medela in Style-Independent home placement test among 76 mothers in USA, August 2012. Medela is a registered trademark of Medela Holding AG
 * BPA Free breast pump: Only associated with the bottle, and other parts that come into contact with breast milk. Following EU regulation, 10/2011