



Philips Avent
Comfort Double electric
breast pump

More comfort, more milk

Breast pump with soft massage cushion

When you are comfortable and relaxed, your milk flows more easily. That is why we created our most comfortable breast pump yet: Sit comfortably with no need to lean forward and let our soft massage cushion gently stimulate your milk flow.



SCF334/16

Sit more comfortably with no need to lean forward

- More comfortable pumping position due to unique design

Gently stimulates let down and milk flow

- Soft massage cushion with warm feel

Other benefits

- Includes our Natural bottle and nipple for natural latch on
- Intuitive assembly. Easy visual matching of parts
- Easy cleaning due to the small number of separate parts

Choose your most comfortable setting

- Features a gentle stimulation mode and 3 pumping settings

More milk in less time

- The perfect time saver for moms

PHILIPS
AVENT

Highlights

More comfortable position



The breast pump has a unique design, so your milk flows directly from your breast into the bottle, even when you are sitting up straight. This means you can sit more comfortably when pumping: no need for you to lean forward to make sure all your milk ends up in the bottle. Sitting comfortably and being relaxed when pumping, naturally helps your milk to flow more easily.

Simple settings to choose from



When switched on, the pump automatically starts in gentle stimulation mode to get your milk flowing. Then choose from 3 pumping settings to make milk flow most comfortable for you.

Easy cleaning

Cleaning is easy, thanks to the small number of separate parts. Your milk will never come into contact with the tubing and base unit. All parts are dishwasher proof, except the electrical parts.

Perfect time saver for moms

Simultaneous pumping at both breasts is proven more efficient and may even boost your ability to produce breast milk.*

Intuitive assembly



Easy visual matching of parts for intuitive assembly.

Includes our Natural nipple



The wide breast shaped nipple promotes natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

Soft massage cushion

Our massage cushion has a soft velvety texture that gives a warm feel to the skin for comfortable, gentle stimulation of your milk flow. During pumping, it gently mimics your baby's suckling for natural let down and clinically proven effectiveness. The soft massage cushion is designed to comfortably fit the majority of moms. Should you need it, a smaller size cushion is available separately.

Specifications

Country of origin

- England

Development stages

- Stage: 0 - 6 months

Material

- Bottle: BPA free*, Polypropylene
- Breast pump: Polypropylene
- Nipple: BPA free*, Silicone

What is included

- Extra soft newborn flow nipple: 4 pcs
- Sealing disc for milk storage: 2 pcs

Design

- Bottle design: Ergonomic shape, Wide neck
- Breast pump design: Compact design

What is included

- Natural bottle 4oz: 4 pcs
- Breast pump body: 2 pcs
- Base unit incl. tubing: 1 pcs
- Large size cushion: 2 pcs
- Storage cups adaptors: 2 pcs
- Breast milk storage cups 6oz: 10 pcs
- Electric breast pump diaphragm: 2 pcs
- Manual breast pump converter: 1 handle and 1 manual breast pump diaphragm
- Breast milk bags: 25 pcs
- Set of 2 comfort breast pad: 1 pcs
- Set of 2 night pads: 1 pcs

Ease of use

- Breast pump use: Intuitive assembly, Easy cleaning, Fully compatible range

Functions

- No leaning forward: Sit in a comfortable position
- Soft massage cushion: Gentle stimulation
- Settings: 3 Expression settings, 1 Stimulation mode
- More milk in less time: The perfect time saver

What is included

- Travel bag: 1 pcs



Issue date 2021-04-08

Version: 5.1.1

© 2021 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* *A randomised controlled trial to compare methods of milk expression after preterm delivery (Jones et al ADC 2001;85:F91).