

# PHILIPS

## AVENT

Comfort Double  
electric breast pump

Natural

Includes 2x 4oz bottle



Product Picture

SCF334/02



## More comfort, more milk

### Breast pump with massage cushion

When you are comfortable and relaxed, your milk flows more easily. That is why we created our most comfortable breast pump yet: Sit comfortably with no need to lean forward and let our soft massage cushion gently stimulate your milk flow.

#### **Sit more comfortably with no need to lean forward**

- More comfortable pumping position due to unique design

#### **Gently stimulates let down and milk flow**

- Soft massage cushion with warm feel

#### **Other benefits**

- Intuitive assembly. Easy visual matching of parts
- Compact lightweight design, includes handy travel bag
- Easy cleaning due to the small number of separate parts
- Easy operation at the touch of a button

#### **Choose your most comfortable setting**

- Features a gentle stimulation mode and 3 pumping settings

#### **More milk in less time**

- The perfect time saver for moms

# Highlights

## More comfortable position



The breast pump has a unique design, so your milk flows directly from your breast into the bottle, even when you are sitting up straight. This means you can sit more comfortably when pumping: no need for you to lean forward to make sure all your milk ends up in the bottle. Sitting comfortably and being relaxed when pumping, naturally helps your milk to flow more easily.

## Simple settings to choose from

When switched on, the pump automatically starts in gentle stimulation mode to get your milk flowing. Then choose from 3 pumping settings to make milk flow most comfortable for you.

## Soft massage cushion



Our massage cushion has a new soft velvety texture that gives a warm feel to the skin for comfortable, gentle stimulation of your milk

flow. The cushion is designed to gently mimic your baby's suckling to help stimulate let down.

## Easy cleaning

Cleaning is easy, thanks to the small number of separate parts. Your milk will never come into contact with the tubing and base unit. All parts are dishwasher proof, except the electrical parts.

## Easy operation

Double electric breast pumps are ideal for moms who pump milk regularly. Save time by comfortably pumping from both breasts at once.

## Perfect time saver for moms

Simultaneous pumping at both breasts is proven more efficient and may even boost your ability to produce breast milk.\*

## Intuitive assembly

Easy visual matching of parts for intuitive assembly.

## Compact design

The breast pump has a compact design, which makes it easy to hold and position on your breast. The small, lightweight base unit can easily be placed within comfortable reach for full control when pumping. For extra convenience in transport and storage, the tube simply wraps around the base unit. The breast pump comes with a handy travel bag.

# Specifications

## Breast pump

**Material:** BPA free\*

## Design

**Breast pump design:** Compact design

**Bottle design:** Ergonomic shape, Wide neck

## Country of origin

England

## Material

**Bottle:** Polypropylene, BPA free\*

**Breast pump:** Polypropylene, BPA free\*

**Nipple:** Silicone, BPA free\*

## What is included

**Extra soft newborn flow nipple:** 2 pcs

**Breast pad sample packs:** 2 (2 Day pads & 2 Night pads) pcs

**Sealing disc for milk storage:** 2 pcs

**Natural bottle 4oz:** 2 pcs

**Standard size cushion:** 2 pcs

**Breast pump body:** 2 pcs

**Travel bag:** 1 pcs

**Base unit incl. tubing:** 1 pcs

**Spare diaphragm:** 2 pcs

**Travel cover:** 2 pcs

## Ease of use

**Breast pump use:** Intuitive assembly, Easy cleaning, Fully compatible range

## Functions

**No leaning forward:** Sit in a comfortable position

**Soft massage cushion:** Gentle stimulation

**Settings:** 3 Expression settings, 1 Stimulation mode

**More milk in less time:** The perfect time saver

## Development stages

**Stage:** 0 - 6 months

