

Philips  
Avent Royal Gift Set

Classic

SCD683



## Good feed, good sleep

### LIMITED EDITION

The royal collection SCD683/31 by Philips Avent including 1 Classic 9 oz feeding bottle and two white classic soothers 6-18m.

#### **Benefits and features:**

- 0% BPA Product

#### **Vents air away from your baby's tummy**

- Clinically proven anti-colic system

#### **More comfortable feed at baby's chosen pace**

- Easy latch on due to the unique valve on the nipple

#### **Reduce fussing and discomfort**

- Clinically proven to significantly reduce fussing

#### **To fit your baby's growing needs**

- Five different nipple flow rates are available

# PHILIPS

# Highlights

## 0% BPA



This product contains 0% BPA

## Unique anti-colic system




As your baby feeds, the unique valve on the nipple flexes to allow air into the bottle instead of your baby's tummy.\*

## Easy latch on

Unique valve on the nipple flexes to your baby's feeding rhythm. Milk will only flow at the pace


chosen by your baby to minimize overeating and spit-up, burping and gas

## Significantly reduces fussing



Sleep and nutrition are vital to your baby's health and happiness. A randomised clinical trial was carried out to see whether infant baby bottle design affects "infant behavior". The Philips Avent Classic baby bottle was shown to significantly reduce fussing by approximately 28 minutes a day as compared to the comparator bottle (46 min vs 74 min,  $p=0.05$ ) This was especially true during the night-time.\*\*

## Fits your baby's growing needs



Five different nipple flow rates are available.

# Specifications

## Development stages

- Stage: 6 months +

## Bottle

- Material: BPA free

## Material

- BPA free\*

## What is included

- Silicone Pacifier: 2 pcs

## Design

- Bottle design: Wide neck, Ergonomic shape

## Features

- Can be sterilized

## Ease of use

- Bottle use: Easy to clean, Easy to hold

## Functions

- Latch on: Easy latch on, Promotes nipple acceptance
- Nipple: Flexes to feeding rhythm

## Functions

- Anti-colic valve: Two piece anti-colic system



Issue date 2021-04-08

Version: 5.0.1

© 2021 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)