

A close-up photograph of a baby's face as they are being fed with a white Philips Avent Natural baby bottle. The baby is wearing a white onesie with blue star patterns. The bottle has the 'AVENT' logo on it. The background is softly blurred, showing a person in a green shirt.

PHILIPS

AVENT

Natural baby bottle

3 Bottles

4oz/125ml

Newborn flow nipple

0m+

 Product Picture

SCF010/37

The most natural way to bottle feed

Natural latch on

Our Natural bottle with an ultra soft nipple more closely resembles the breast. The wide breast-shaped nipple with flexible spiral design and comfort petals allows natural latch on and makes it easy to combine breast and bottle feeding.

Easy to combine breast and bottle feeding

- Soft nipple designed to mimic the feel of the breast
- Natural latch on due to the wide breast shaped nipple

Vents air away from your baby's tummy

- Unique anti-colic valve technology

Other benefits

- This bottle is BPA free*
- Simple to use and clean, quick and easy assembly
- Ergonomic shape for maximum comfort
- Compatible with the Philips Avent Natural range
- Nipples with different flow rates available

A comfortable and contented feed for your baby

- Flexible spiral design combined with comfort petals

Highlights

Natural latch on



The wide breast shaped nipple promotes a natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

Anti-colic valve technology



Our unique anti-colic valve technology was designed to reduce colic and discomfort, by venting air away from your baby's tummy.

Simple to use and clean




Wide bottle neck makes filling and cleaning easy. Only a few parts for quick and simple assembly.

Ultra soft nipple



The nipple has an ultra soft texture, designed to mimic the feel of the breast.

Choosing the right nipple



The Philips Avent Natural range offers different nipple softness and increasing flow rates for every development stage of your baby.

Compatibility across the range



Mix and match our breast pump, bottle and cup parts, and create the product that works for you, when you need it!

Flexible spiral design



Flexible spiral design, combined with our unique comfort petals to create a flexible nipple, allowing for a more natural feed without nipple collapse.

Ergonomic shape



Due to the unique shape, the bottle is easy to hold and grip in any direction for maximum comfort, even for baby's tiny hands.

This bottle is BPA free*



The Philips Avent Natural bottle is made of BPA free* material (polypropylene).

Specifications

Design

Bottle design: Ergonomic shape, Wide neck

Material

Bottle: Polypropylene, BPA free*

Nipple: BPA free*, Silicone

What is included

Baby bottle: 3 pcs

Ease of use

Bottle use: Easy to hold, Easy to assemble, Easy to clean

Functions

Anti-colic valve: Advanced anti-colic system

Latch on: Natural latch on, Easy combine breast and bottle

Nipple: Unique comfort petals, Ultra soft and flexible nipple

Development stages

Stage: 0-6 months



* 0% BPA, following EU regulation 10/2011

* What colic is, and how it affects babies? Colic is caused in part by swallowing air while feeding, which creates discomfort in a baby's digestive system. Symptoms include crying, fussing, gassiness and spit-up.