



Philips Avent  
Double electric breast  
pump

Express more milk in less time  
Includes soft massage cushion  
Travel bag and 2 bottles



SCF334/31

## More comfort, more milk, in less time\*

### Double electric breast pump with massage cushion

The Ultra Comfort Double Breast Pump is ideal if you don't have much time. Reduce your expression time & sit comfortably as the soft massage cushion gently stimulates your flow. The pump is quiet & easy to set up, personalise, use & clean.

#### Easy to set up, express and clean

- Few separate parts and intuitive design
- Closed system, designed for hygienic expressing
- Includes our Natural bottle and nipple for natural latch on

#### Quiet pump for discreet expression, anytime

- Designed to let you express quietly

#### 1 stimulation mode and 3 expression settings

- Choose the setting that is effective & comfortable for you

#### Gently stimulates milk flow

- Soft massage cushion with massaging petals
- More milk in less time with double pumping

#### No need to lean forward

- Comfortable pumping position due to unique design

**PHILIPS**  
**AVENT**

Double electric breast pump  
Express more milk in less time Includes soft massage cushion, Travel bag and 2 bottles

SCF334/31

## Highlights

### Relaxing expression position



The breast pump has a unique design, so your milk flows directly from your breast into the bottle or container, even when you are sitting up straight. This means you can sit comfortably when pumping: no need for you to lean forward to make sure all your milk ends up in the bottle. Sitting comfortably and being relaxed when pumping, helps your milk to flow more easily.

### Simple settings to choose from

When switched on, the pump automatically starts in gentle stimulation mode to get your milk flowing. Then choose from 3 pumping settings to make milk flow comfortable for you.

### Soft massage cushion



Our massage cushion has a soft velvety texture that gives a warm feel to the skin for

comfortable, gentle stimulation of your milk flow. The iconic petal cushion is designed to mimic baby suckling, intended to gently help stimulate let down.

### More milk in less time



Reduce expression time with double pumping compared to single pumping.

### No milk spillage in the tubes



No milk spillage in the tubes.

### Easy set up, express and clean



Few separate parts and intuitive design. Easily express milk for your baby, anytime.

### Includes our Natural bottle



The wide breast shaped nipple promotes natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

### Quiet pump

Designed for discreet expression, anytime.

Double electric breast pump  
Express more milk in less time Includes soft massage cushion, Travel bag and 2 bottles

SCF334/31

# Specifications

## Country of origin

- England

## Development stages

- Stage: 0 - 6 months

## Material

- Bottle: BPA free\*, Polypropylene
- Nipple: BPA free\*, Silicone
- Breast pump: BPA free\* (food contact parts only)

## What is included

- Sealing disc for milk storage: 2 pcs

## Design

- Bottle design: Ergonomic shape, Wide neck
- Breast pump design: Compact design

## What is included

- Natural bottle 4oz: 2 pcs

- Breast pump body: 2 pcs
- Base unit incl. tubing: 1 pcs
- Travel cover: 2 pcs
- Breast pad sample packs: 1 pack (2pcs)

## Ease of use

- Breast pump use: Easy cleaning, Fully compatible range, Intuitive assembly

## Functions

- No leaning forward: Sit in a comfortable position
- Soft massage cushion: Gentle stimulation
- Settings: 1 Stimulation mode, 3 Expression settings

## What is included

- Travel bag: 1 pcs
- Standard size cushion(19.5 mm): 2 pcs (a larger size cushion is available and sold separately)



Issue date 2021-10-13

Version: 6.1.1

EAN: 08 71010 39535 17

© 2021 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)

\* 83% of breastfeeding mothers agree that Philips Avent Comfort double electric breast pump is comfortable during usage – Independent home placement test among 81 mothers in the USA, August 2012.  
\* More milk: Independent research has shown that there may be a link between higher comfort levels and increased milk production. See [www.philips.com/AVENT](http://www.philips.com/AVENT)  
\* BPA Free breast pump: Only associated with the bottle, and other parts that come into contact with breast milk. Following EU regulation, 10/2011