

PHILIPS

Silicone Frame

Wisp

Replacement part



HH1001

Clear unobtrusive silicone frame

The clarity of the silicone frame draws less attention to the face.

Silicone frame

- Clear silicone frame

Highlights

Unobtrusive design

The clarity of the silicone frame draws less attention to the face.

Specifications

Service

Warranty: 90-days

Design

Compact: Easy fitting and adjustment

Comfort: Elastic crown strap

Patented design: Tip-of the nose cushion design

Quiet exhalation: 360 degrees rotating elbow

Silicone frame: Clear unobtrusive frame

Specifications

Frame material: Silicone/Polycarbonate

