

PHILIPS

Gel Pillows Cushion
(S)

DreamWear

Small size



HH1125



Like wearing nothing at all*

With the freedom and comfort of DreamWear Gel Pillows, you can sleep how you want.

Easy to use

- Features just for you

Innovative gel pillows design

- Easy to fit
- Comfort and performance
- Modular design

Highlights

Easy to fit



Choose the cushion so the pillows cushion blue gel bases sit comfortably under the nose without any gaps. Only the soft silicone tip should be inside your nostrils. Do not push any part of the blue gel base inside your nose.

Innovative design



Less discomfort on the face, nose and nostrils. Users rated DreamWear Gel Pillows more comfortable than their prescribed mask.

Innovative design



Easy to use with one modular frame for all cushion sizes.

Modular design



Modular design for quick assembly and cleaning, with interchangeable pillows and nasal cushions.

Specifications

Product details

Materials: Polycarbonate, Polyurethane gel, Silicone

Warranty: 90 days

Cleaning: Cushion

Hand wash with mild detergent: Rinse thoroughly, air dry¹

Replacement

Inspect daily for wear: Replace at first sign of wear

