

PHILIPS


InfraCare infrared
lamp

650 W

Half body treatment

63.5 cm high

HP3641

 Product Picture

Effective pain relief

Deep penetrating infrared warmth

Philips InfraCare gives effective pain relief for muscles and joints. Its 650 W infrared warmth is very comfortable and penetrates deep into the skin, stimulates blood circulation and gives effective relief for half-body areas of 60x40 cm.

Easy to use

- Adjustable angle with rotation

Effective pain relief

- Deep penetrating Infrared warmth

Half-body treatment

- 650 Watt infrared halogen lamp

Highlights

650 Watt infrared halogen lamp



The Philips InfraCare HP 3643 has been optimised for effective relief of pains in areas covering 60 x 40 cm, such as the entire back, both shoulders and neck or thigh. The InfraCare features innovative infrared halogen lamp technology. The special optics, filter and the powerful halogen lamp of 650 Watt have been optimised, enabling half-body treatment (60 x 40 cm), and much more comfortable warmth, as the infrared warmth is more equally divided over the treatment area.

Adjustable angle with rotation

The InfraCare HP 3643 is easy to use for half-body treatment and for different situations (e.g. lying on a bed or sofa) due to its maximum flexibility. In addition to the height adjustment,

the lamp house can easily be positioned for optimum treatment of the right body part by moving it into a horizontal or vertical position (left/right 90 degrees) and by tilting it up or down (45 degrees).

Deep penetrating Infrared warmth



Infrared light has been proven to effectively relieve aches and pains caused by muscular pains and stiff joints. The comfortable warmth created by infrared light penetrates deep into your skin, stimulates blood circulation and warms your muscles. As your muscles are soothed, they automatically loosen up and relax. Because it makes the tissue more pliable, warmth also reduces stiffness and makes joints more flexible.

Specifications

Easy to use

2 handgrips: for easy transport and positioning

On/off switch

Extendable height: 45-80 cm

Easily positionable lamphouse (left/right): 90

Easily positionable lamphouse (up/down): 45

Half-body treatment

Infrared halogen lamp: 650 W

Treatment area: 60 x 40 cm

Safety

IEC certified: Complies to IEC 60601 and 60335

Effective pain relief

Deep penetrating warmth

Treatment control

Digital timer: with auto shut-off

Dimmer with 3 warmth settings

Pause function

Technical specifications

Voltage: 230 V

Cord length: 300 m

Power: 650 W

Frequency: 50/60 Hz

Insulation: Class II (double isolation)

Lifetime of lamps: 500 hour(s)

Weight and dimensions

Product dimensions: 29.0 x 63.5 x 29.0 (WxHxD) cm

Product weight: 5.5 kg

Qnt. on Euro pallet: 18 pcs

F-box dimensions (= A-box): 33.9 x 80.0 x 36.5 (WxHxD) cm

F-box weight (= A-box): 6.5 kg

Logistic data

CTV code: 884364101000

Country of origin: Germany

Easy to use

Digital user interface: for feedback on time and warmth intensity

Medical appliance

Medical Device Directive: 2007/47/EC, MDD 93/42/EEC

