

PHILIPS

AVENT

Manual breast pump with 3 cups

Easy for on-the-go use

Pump with soft massage cushion

Includes 3 storage cups

Extra soft newborn flow nipple

SCF330/13



More comfort, more milk, easy on the go*

Manual breast pump with massage cushion

Sit comfortably without leaning forward as our soft massage cushion gently stimulates your milk flow. The Avent manual breast pump has few parts, is easy to assemble, use & clean. Lightweight and compact, it's easy to express on-the-go.

Sit more comfortably with no need to lean forward

- More comfortable pumping position due to unique design

Gently stimulates let down and milk flow

- Soft massage cushion with warm feel

Other benefits

- Intuitive assembly. Easy visual matching of parts
- Easy cleaning thanks to the small number of separate parts
- Easy manual operation
- Compatible with other Philips Avent feeding products

Express, store and feed in one

- Includes 3 versatile milk storage cups

Easy to use, store and transport

- Compact lightweight design

Highlights

More comfortable position



The breast pump has a unique design, so your milk flows directly from your breast into the bottle, even when you are sitting up straight. This means you can sit more comfortably when pumping: no need for you to lean forward to make sure all your milk ends up in the bottle. Sitting comfortably and being relaxed when pumping, naturally helps your milk to flow more easily.

Soft massage cushion



Our massage cushion has a new soft velvety texture that gives a warm feel to the skin for comfortable, gentle stimulation of your milk flow. The cushion is designed to gently mimic your baby's suckling to help stimulate let down.

Express, store and feed in one



Use the same cup to express, store and feed breast milk to your baby. Ideal for storage in fridge or freezer

Compact design



The breast pump has a compact design, which makes it easy to position. Its ergonomically shaped handle offers a comfortable hold and full control when expressing milk. Being small and lightweight means it is easy to store and transport, making pumping on-the-go more discreet.

Intuitive assembly



Easy visual matching of parts for intuitive assembly.

Fully compatible range



This breast pump can be used in combination with other feeding products in the Philips Avent range, including our Classic bottles and milk storage containers. Philips Avent also provides a range of breast care accessories to help you breastfeed for longer and enhance your comfort.

Easy cleaning

Cleaning is easy, thanks to the small number of separate parts. All parts are dishwasher proof.

Easy manual operation

The manual breast pump is ideal if you express milk occasionally and value compactness. Easily operated with one hand.

Specifications

Design

Breast pump design: Compact design

Bottle design: Ergonomic shape, Wide neck

Material

Nipple: BPA free*, Silicone

Milk storage cup: BPA free*, Polypropylene

Breast pump: BPA free* (food contact parts only)

What is included

Extra soft newborn flow nipple: 1 pcs

Breast pump body with handle: 1 pcs

Milk storage cup adaptor: 1 pcs

Standard size cushion: 1 pcs

Breast pad sample packs: 1 pack (2pcs)

Milk storage cups: 3 pcs

Ease of use

Breast pump use: Intuitive assembly, Easy cleaning, Fully compatible range

Functions

No leaning forward: Sit in a comfortable position

Soft massage cushion: Gentle stimulation

Includes milk storage cups: Express, store and feed in one

Development stages

Stage: 0 – 6 months

