

PHILIPS

IntraSight

Interventional
applications platform

Mobile

Smart. Simple. Seamless.



Smart

IntraSight is an interventional platform protected by Windows 10, your best defense against cybersecurity threats.



Simple

Versatility from a touchscreen, modern and easy to maneuver mobile cart design.



Seamless

Full bedside procedural control with TSM in sterile field.

Uniting today's best-in-class imaging on a secure applications based platform ready for tomorrow's innovations.



TSM tableside control


Philips allows you to run an entire case without breaking scrub. With our intuitive interface and new touch screen module (TSM), you get total control within the sterile field. Run an entire case tableside as you quickly navigate, annotate and measure to your exact specifications. All with the ease and responsiveness of using a tablet, even when wearing gloves and under a sterile drape.



Configuration

Series 3

Interventional applications platform	•
Imaging (IVUS) license	•
Includes ChromaFlo (IVUS)	Available on Visions PV .014 and .018 and Reconnaissance IVUS catheters only
Touch Screen Module (TSM)	Optional
Philips Remote Services	•
Demonstration mode to facilitate training	•



Simplified annotation by procedure type

Delivering an outstanding user experience with enhanced annotation capabilities specific to venous, arterial, AV access, thoracic and abdominal procedures.

