

Philips Philips Avent  
Breast pump motor unit  
and handle

**Breast pumps**  
Single electric breast pump

SCF932/01



## Pumping and suction adjustment

Check for compatibility below

The breast pump motor unit learns your desired rhythm after two presses and then repeats the rhythm automatically. The handle allows you to adjust the pressure, so you can find the suction that's right for you.

### **Ease of using Philips spare parts**

- Easily renew your product with original Philips parts

**PHILIPS**

Breast pump motor unit and handle  
Breast pumps Single electric breast pump

SCF932/01

# Specifications

# Highlights

## Replaceable parts

- Fit product types: SCF292, SCF302, SCF312

## Easy replacement of product parts

From time to time your product needs a facelift, and with Philips consumer replacement parts to renew your product, it's never been so easy! All this with guaranteed Philips quality.



Issue date 2024-04-13

Version: 4.4.1

© 2024 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)