

PHILIPS



Full Face Headgear

DreamWear

Full Face Replacement Part
Standard Headgear



HH1131

Like wearing nothing at all.*

With the freedom and comfort of DreamWear Full Face, you can sleep how you want.

Magnetic clips

- The headgear clips and mask frame contain magnets

Comfort and ease of use

- Soft fabric headgear with magnet clips

Highlights

Soft and comfortable



The soft fabric headgear is comfortable to wear at night while the magnet clips make it easy to adjust for the right fit.

Magnetic clips



Some medical implants may be contraindicated for use with magnetic devices and this mask may not be suitable for use. Contact your healthcare professional before you use this mask.

Specifications

Warning: Magnetic interference

Beware with pacemakers: Read IFU's before use.

Replacement

Inspect daily for wear: Replace at first sign of wear

Product details

Materials: Nylon, Polyurethane foam, Spandex

Warranty: 90 days

