

PHILIPS

Frame DreamWear

Replacement Part

Standard frame

Available in multiple sizes



HH1117/00



Like wearing nothing at all.*

With the freedom and comfort of DreamWear, you can sleep how you want.

Easy to use


- Features just for you

Innovative design

- Soft, flexible frame
- Air tube connection at the top of the head

Highlights

Easy to use frame



Easy to use, choose one of the three frame size that fits you best. All the frame sizes fit all the cushion sizes.

Innovative frame design



The soft, flexible frame is comfortable to wear at night. Users have stated it feels like wearing nothing at all.*

Freedom of movement



By moving the tube location to the top of the head, DreamWear offers more freedom of movement during the night. Users stated they didn't need to choose a sleep position when wearing DreamWear.*

Specifications

Innovative design

Comfort and performance: Soft fabric straps

Product details

Materials: Polycarbonate, Silicone

Warranty: 90 days

Cleaning: Frame

Hand wash with mild detergent: Rinse thoroughly, air dry¹

Replacement

Inspect daily for wear: Replace at first sign of wear



* * Data analysis of Aug 2017 Patient Preference trial where n = 85 and prescribed masks include Resmed Airfit F10 and F20, Respiroics Amara View and F&P Simplus; based on 10 days of use.

* *2017 Philips User Preference Questionnaire

* ¹Please refer to Instructions For Use for more detailed cleaning instructions.