

# PHILIPS

Full Face Headgear  
(Standard)

DreamWear

Full Face Replacement Part  
Standard Headgear



HH1131



## Like wearing nothing at all.\*

With the freedom and comfort of DreamWear Full Face, you can sleep how you want.

### **Magnetic clips**

- The headgear clips and mask cushion contain magnets

### **Comfort and ease of use**

- Soft fabric headgear with magnet clips

# Highlights

## Soft and comfortable



The soft fabric headgear is comfortable to wear at night while the magnet clips make it easy to adjust for the right fit.

## Magnetic clips



Contact your healthcare professional before you use this mask. Some medical devices may be affected by magnetic fields. The magnetic clips in this mask should be kept at least 2 in. (50mm) away from any active medical device with special attention to implanted devices such as pacemakers, defibrillators and cochlear implants. Read in-pack IFU's for more information.

# Specifications

## Warning: Magnetic interference

Beware with pacemakers: Read IFU's before use.

## Replacement

Inspect daily for wear: Replace at first sign of wear

## Product details

Materials: Nylon, Polyurethane foam, Spandex

Warranty: 90 days

