

# PHILIPS

Under-the-nose  
cushion

DreamWear Full  
Face

Comfortably sits under the nose

Easily clips to headgear with  
magnets

Prevents red marks

Available in multiple sizes



HH1128/02



## Like wearing nothing at all.\*

With the freedom and comfort of DreamWear Full Face, you can sleep how you want.

#### **Innovative design**

- A modular design for more options

#### **Comfort and ease of use**

- Easy to use
- Under the nose full face cushion

# Highlights

## Easy to use



Only a few parts for easy assembly, disassembly and cleaning.

## Freedom of choice



The modular design allows users to choose the cushion that works best for them, all on one frame.

## No more nose bridge red marks



The under the nose innovative design of the cushion means no more red marks, discomfort, and irritation on the nose bridge.\*\*

# Specifications

## Product details

**Materials:** Polycarbonate, Silicone

**Warranty:** 90 days

## Cleaning: Cushion

**Hand wash with mild detergent:** Rinse thoroughly, air dry<sup>1</sup>

## Replacement

**Inspect daily for wear:** Replace at first sign of wear



\*\* Data analysis of Aug 2017 Patient Preference trial where n = 85 and prescribed masks include Resmed Airfit F10 and F20, Respiration Amara View and F&P Simplus; based on 10 days of use.

<sup>1</sup> Please refer to Instructions For Use for more detailed cleaning instructions.