

PHILIPS

Under-the-nose
cushion

DreamWear Full
Face

Comfortably sits under the nose

Easily clips to headgear with
magnets

Prevents red marks

Available in multiple sizes



HH1128/02



Like wearing nothing at all.*

With the freedom and comfort of DreamWear Full Face, you can sleep how you want.

Innovative design

- A modular design for more options

Comfort and ease of use

- Easy to use
- Under the nose full face cushion

Highlights

Easy to use



Only a few parts for easy assembly, disassembly and cleaning.

Freedom of choice



The modular design allows users to choose the cushion that works best for them, all on one frame.

No more nose bridge red marks



The under the nose innovative design of the cushion means no more red marks, discomfort, and irritation on the nose bridge.**

Specifications

Cleaning: Cushion

Hand wash with mild detergent: Rinse thoroughly, air dry¹

Replacement

Inspect daily for wear: Replace at first sign of wear

Product details

Materials: Polycarbonate, Silicone

Warranty: 90 days



* * Data analysis of Aug 2017 Patient Preference trial where n = 85 and prescribed masks include Resmed Airfit F10 and F20, Respironics Amara View and F&P Simplus; based on 10 days of use.

¹Please refer to Instructions For Use for more detailed cleaning instructions.