

PHILIPS

Full Face Mask - Fit
Pack (S, M, L)

Amara View

Fit Pack - 3 cushion sizes (S, M,
L)

1 modular frame



HH1231



Enjoy the view

With its innovative under-the-nose design, Amara View offers an open field of vision and prevents red marks, discomfort or irritation on the nose bridge.

Open field of vision

- Widest field of vision of all leading full face masks

Innovative design

- Headgear with crown strap
- Four-part modular design
- Innovative design

The smallest and lightest choice

- The smallest* and lightest** of all leading full face masks

Highlights

Smallest and lightest choice



Smaller and lighter than traditional full face masks, Amara View covers less of the patient's face than comparable masks and eliminates the bulky cushion and frame in front of the eyes.

Unobstructed view



It's easy to wear glasses, read, watch TV, and use a computer or tablet before falling asleep.

Headgear with crown strap



Headgear with soft and comfortable fabric straps. The crown strap provides quick and easy fitting adjustment.

Four-part modular design



Modular design for quick assembly and cleaning, with interchangeable cushions available in four sizes.

Under-the-nose cushion



Prevents red marks, discomfort or irritation on the nose bridge

Specifications

Innovative design

Wide field of vision: Wear glasses, read or watch TV

Specifications

Operating pressure: 5 -30 cmH₂O

Sound pressure level: 24.2 dBA

No. of parts: Four

Cushion material: Silicone

Headgear material: UBL/mesh blend

Cleaning instructions: mask: Rinse thoroughly, air dry

Cleaning: headgear, tubing: Hand wash with mild detergent, Rinse thoroughly, hang to dry

Headgear clips/Tubing material:

Thermoplastic polyester elastomer

Quick-release feature material:

Polyoxymethylene (Acetal)

Frame/Elbow/Swivel material:

Polycarbonate

Intent of use

To be used on patients <gt;/>30 kg: with CPAP prescription

Contraindications:

Impaired cough reflex: impaired cardiac sphincter

Recent eye surgery, dry eyes,: hiatal hernia, major reflux

Service

Warranty: 90-days

