

PHILIPS

Mask

Nuance

Standard

Pack

3 pillows (S, M, L)



HH1100



A new choice for better sleep

The nasal pillow technology of Nuance delivers exceptional comfort and seal, helping you to support better sleep.

Easy to use

- Light-weight, flexible tubing for less pull on the mask.

Appealing, minimalistic design

- Simple design is easier to use and fit.

Gel pillows

- The pillows fit difference nostril sizes.

Greater stability

- Stays in place throughout the night.

Highlights

Lightweight and simple



Easier to assemble, disassemble and use than the leading pillows mask.

Provides comfort and seal



Nuance gel pillows cushion is designed differently to provide a more comfortable experience. The gel base of the pillows should fit comfortably under your nose.

Innovative non-slip headgear



Ideal for patients that want to sleep freely.

Patient experience



During our patient preference study the Nuance platform was preferred over the leading nasal pillows mask for seal, comfort, ease of use and for less nasal irritation.

Specifications

Specifications

Operating pressure: 4 -20 cmH2O

Sound pressure level: 19dBA

No. of parts: Four

Cushion material: Silicone rubber

Headgear material: Nylon UBL/Urethane foam/Lycra

Frame material: Nylon, Lycra, Urethane foam and Trico mesh

Cleaning instructions: mask: Hand wash with mild detergent, Rinse thoroughly, air dry

Cleaning: headgear, tubing: Hand wash with mild detergent

Swivel/Tubing/Elbow material:

Polypropylene

Pillows material: Urethane gel

Intent of use

To be used on patients <gt/>30 kg: with CPAP prescription

Design

Easy to use: Easy assembly and disassembly

Comfort: Fabric frame

Service

Warranty: 90-days

