

**PHILIPS**

**Mask**


**Nuance**

Standard

Pack

3 pillows (S, M, L)

HH1100



# Highlights

## Lightweight and simple



Easier to assemble, disassemble and use than the leading pillows mask.

## Provides comfort and seal



Nuance gel pillows cushion is designed differently to provide a more comfortable experience. The gel base of the pillows should fit comfortably under your nose.

## Innovative non-slip headgear



Ideal for patients that want to sleep freely.

## Patient experience



During our patient preference study the Nuance platform was preferred over the leading nasal pillows mask for seal, comfort, ease of use and for less nasal irritation.

# Specifications

## Service

Warranty: 90-days

## Design

Easy to use: Easy assembly and disassembly

## Design

Comfort: Fabric frame

## Specifications

Operating pressure: 4 -20 cmH2O

Sound pressure level: 19dBA

N° of parts: Four

Cushion material: Silicone rubber

Headgear material: Nylon UBL/Urethane foam/Lycra

Frame material: Nylon, Lycra, Urethane foam and Trico mesh

Cleaning instructions: mask: Hand wash with mild detergent, Rinse thoroughly, air dry

Cleaning: headgear, tubing: Hand wash with mild detergent

Swivel/Tubing/Elbow material:

Polypropylene

Pillows material: Urethane gel

## Intent of use

To be used on patients <gt/>30 kg: with CPAP prescription

