

# PHILIPS

## AVENT

### Manual breast pump with bottle

Easy for on-the-go use

Includes 4oz / 125ml bottle

Includes bottle and nipple

Travel cover and sealing disc



SCF330/20



## More comfort, more milk, easy to take on-the-go\* Manual breast pump with massage cushion

Sit comfortably without leaning forward as our soft massage cushion gently stimulates your milk flow. The Avent manual breast pump has few parts, is easy to assemble, use & clean. Lightweight and compact, it's easy to express on-the-go.

#### **Easy to combine breast and bottle feeding**

- Includes our Natural bottle and nipple for natural latch on

#### **Easy to set up, express and clean**

- Few separate parts and intuitive design
- Compact lightweight design

#### **Clinically proven less leaning forward\***

- Relaxing position, thanks to its special design

#### **Soft massage cushion with massaging petals**

- Soft massage cushion gently stimulates milk flow

# Highlights

## Relaxing position



The breast pump has a unique design, so your milk flows directly from your breast into the bottle, even when you are sitting up straight. This means you can sit more comfortably when pumping: no need for you to lean forward to make sure all your milk ends up in the bottle. Sitting comfortably and being relaxed when pumping, helps your milk to flow more easily.

## Soft massage cushion



Our massage cushion has a soft velvety texture that gives a warm feel to the skin for comfortable, gentle stimulation of your milk flow. The iconic petal cushion is designed to mimic baby suckling, intended to gently help stimulate let down.

## Easy on-the-go



The breast pump is small and lightweight, which means it is easy to store and transport, making pumping on-the-go more discreet.

## Easy set up, express and clean



Few separate parts and intuitive design. Easily express milk for your baby, anytime.

## Includes our Natural bottle



The wide breast shaped nipple promotes natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

# Specifications

## Design

**Breast pump design:** Compact design  
**Bottle design:** Ergonomic shape, Wide neck

## Material

**Bottle:** BPA free\*, Polypropylene  
**Nipple:** BPA free\*, Silicone  
**Breast pump:** BPA free\* (food contact parts only)

## What is included

**Extra soft newborn flow nipple:** 1 pcs  
**Breast pump body with handle:** 1 pcs  
**Sealing disc for milk storage:** 1 pcs  
**Natural bottle 4oz:** 1 pcs  
**Travel cover:** 1 pcs  
**Breast pad sample packs:** 1 pack (2pcs)  
**Standard massage cushion:** 1 pcs(19.5mm, a larger size cushion is available separately)

## Ease of use

**Breast pump use:** Easy cleaning, Fully compatible range, Intuitive assembly

## Functions

**No leaning forward:** Sit in a comfortable position  
**Soft massage cushion:** Gentle stimulation

## Development stages

**Stage:** 0 - 6 months

