



A collection of AVEEIT breastfeeding products. In the foreground, there is a manual breast pump with a large white dome and a purple base. Behind it is a double electric breast pump with a purple top and two white bottles. To the right, there are several white bottles of different sizes, some with blue caps, and a small tube of AVEEIT cream. In the background, there is a purple box with the AVEEIT logo. In the foreground, there are two white breast nipples and a purple ribbon.

Breastfeeding made easy

The Comfort Breastfeeding Support Kit is the ideal set for mums who plan to express their milk. The set contains all you need to pump, store and feed your breast milk.

Express - Single Electric Breast Pump

- Choose the most effective and comfortable setting for you
- Designed to let you express quietly
- Soft massage cushion gently stimulates milk flow
- Few separate parts and intuitive design
- Power cushion for extra suction
- Comfortable pumping position due to unique design

Feed - Natural bottle and teat

- Includes our Natural bottle and teat for natural latch-on

Store - Storage cups

- For safe storage and transport
- Pump in storage containers

Care - Thermogel and disposable breast pads

- Silky soft feel, breathable material

Highlights

Relaxing expression position



The breast pump has a unique design, so your milk flows directly from your breast into the bottle or container, even when you are sitting up straight. This means you can sit comfortably when pumping: there's no need for you to lean forward to make sure all your milk ends up in the bottle. Sitting comfortably and being relaxed when pumping helps your milk to flow more easily.

Soft massage cushion



Our massage cushion has a soft, velvety texture that gives a warm feel to the skin for comfortable, gentle stimulation of your milk flow. The iconic petal cushion is designed to mimic baby suckling, intended to gently help stimulate let-down.

Soft power cushion



Choose the power cushion for the moments you need a stronger suction.

Simple settings to choose from

When switched on, the pump automatically starts in gentle stimulation mode to get your milk flowing. Then choose from 3 pumping settings to make milk flow in the way that is most comfortable for you.

Quiet pump

Designed for discreet expression.

Easy to set up, express and clean



Few separate parts and intuitive design. Easily express milk for your baby, anytime.

Includes our Natural bottle



The wide, breast-shaped teat promotes natural latch-on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

Pump in storage containers



With the milk storage cup adapter, you can easily pump in the milk storage container. No leakage, no spillage of your precious breast milk.

Leak-proof lid



Philips Avent storage cups have a lid that creates a secure seal for safe storage and transport.

Silky-soft feel



The Philips Avent disposable breast pads are made of silky soft topsheet and breathable material. They have multiple layers for super absorbency.

Specifications

Design

Bottle design: Ergonomic shape, Wide neck

Breast pump design: Compact design

Country of origin

England

Material

Bottle: BPA-free*, Polypropylene

Teat: BPA-free*, Silicone

Breast pump: BPA-free* (food contact parts only)

Breast pads: Natural breathable material, Silky soft top-sheet

Storage cups: BPA free*, Polypropylene

What is included

Milk storage cup adapter: 2 pcs

Sealing disc for milk storage: 1 pcs

Natural bottle 125 ml/4 oz: 1 pcs

Natural bottle 9 oz: 1 pcs

Breast pump body: 1 pcs

Base unit incl. tubing: 1 pcs

Milk storage cup (180 ml/6 oz): 10 pcs

Washable pouch: 1 pcs

Teat cream: 1 pcs

Disposable breast pads: 18 pcs

Massage cushion (19.5 mm): 1 pcs (a larger size cushion is available and sold separately)

Power cushion (27mm): 1 pcs

Teat protectors: 2 pcs

Ribbon bra clip: 1 pcs

Ease of use

Breast pump use: Easy cleaning, Fully compatible range, Intuitive assembly

Functions

No leaning forward: Sit in a comfortable position

Soft massage cushion: Gentle stimulation

Settings: 3 Expression settings, 1 Stimulation mode

Development stages

Stage: 0–6 months

