

# PHILIPS

## AVENT

### Manual Breast Pump

Easily pump and feed

Standard set



SCF430/10



## Comfortable pumping, wherever you go

Natural Motion Technology, for quick milk flow

Enjoy the gentle comfort of the Philips Avent portable manual breast pump. With Natural Motion Technology, inspired by baby's own suckling action for a quick letdown. You can easily adjust rhythm and vacuum. Fits practically all nipples.

#### Easy to clean

- Easy to clean and set up

#### Comfortable expression

- Express and feed
- Easily triggers milk release before or between feeding
- Express without leaning forward
- Choose your rhythm, follow your milk flow

#### Portable and convenient

- Ideal for moms who pump on-the-go

#### Simple, gentle and effective

- One-size soft adapting silicone cushion
- Helps you release milk fast

# Highlights

## Natural Motion Technology



The innovative silicone flange – inspired by baby’s unique suckling action – stimulates the nipple to help start your milk flow, fast. For comfortable and effective expression.

## Adapts to your nipples



One size fits all. Because we all come in different shapes and sizes, the silicone cushion gently flexes and adapts to fit your unique anatomy. Fits 99.98% of all nipple sizes\* (up to 30mm).

## Portable and easy to use



The breast pump is small and lightweight, which means it is easy to store and transport, making pumping on-the-go simple and discreet.

## Fast milk release




A unique, portable manual breast pump – ideal for mothers who need to gently release their milk while away from home, or before breastfeeding.

## Suction and flow control



Feel free to adjust the pumping rhythm and suction to suit your needs, comfort and milk flow.

## Sit comfortably



Relax and sit comfortably while you express – thanks to our design that lets you sit upright, rather than having to lean forward.

## Easy to clean and set up



Because there are only a few parts, our new pump is easy to set up and clean.

## Express and feed



To feed using your expressed milk, simply attach the soft nipple to the bottle container. To store your milk for feeding later on, close the bottle using the sealing disc.

# Specifications

## Material

**Bottle:** Polypropylene, BPA free\*

**Teat:** Silicone, BPA free\*

**Nipple:** BPA free

**Breast pump:** BPA free\*

## What is included

**4oz Bottle with Om+ nipple:** 1 pcs

**Disposable breast pads:** 2 pcs

**Sealing disc:** 1 pcs

**Cushion cover:** 1 pcs

**Expression kit with handle:** 1 pcs

## Functions

**No leaning forward:** Sit in a comfortable position

**Soft & adaptive cushion:** Gentle stimulation

**Portable and lightweight:** Easy on the go



\* Based on: (1)Mangel et al. Breastfeeding difficulties, breastfeeding duration, maternal body mass index, and breast anatomy: are they related?. Breastfeeding Medicine, 26 April 2019, ((109 participants, Israel); (2)Ziemer et al. Skin changes and pain in the nipple during the 1st week of lactation.

\* Journal of Obstetric, Gynecologic & Neonatal Nursing, May 1993, (20 participants (Caucasian), USA); (3)Ramsay et al. Anatomy of the lactating human breast redefined with ultrasound imaging, 2005, (28 participants, Australia).

\* BPA Free breast pump: Only associated with the bottle, and other parts that come into contact with breast milk. Following EU regulation, 10/2011