

PHILIPS

Efficia CM

Patient monitor



Efficia patient monitor product structure



Model
Part number
Screen

CM10
863301
10"

CM12
863303
12"

CM100
863300
10"

CM120
863302
12"

CM150
863304
15"

Measurement parameter options (accessories must be ordered separately)

3- and 5-lead ECG		•	•	•	•	•
QT/QTc Interval		•	•	•	•	•
12-lead ECG	B07		□		□	□
Philips ST analysis		•	•	•	•	•
Impedance respiration		•	•	•	•	•
Philips SpO ₂ ¹	B01	□	□	□	□	□
Masimo rainbow SET ²	B02	□	□	□	□	□
Nellcor OxiMax SpO ₂ ²	B11	□	□	□	□	□
Masimo rainbow SpHb/SpOC ³	B03			□	□	□
Masimo rainbow SpCO ³	B04			□	□	□
Masimo rainbow RRA ³	B05			□	□	□
Philips NBP		•	•	•	•	•
Dual continuous temp		•	•	•	•	•
Dual IBP with overlapp IBP waveforms and Pulse pressure variation - PPV	B10	□	□	□	□	□
CO ₂ Oridion Microstream ⁴	B20	□	□	□	□	□
CO ₂ Respiration Capnostat Mainstream ⁴ (Efficia CM patient monitor configured with the B21 option are compatible with the Respiration Loflo sidestream sensor (M2741A). M2741A must be ordered separately.)	B21	□	□	□	□	□
CO ₂ Respiration CapnoTrak Sidestream ⁴	B22	□	□	□	□	□
Cardiac output ⁵ (with hemodynamics calculations and wedge pressure)	B30		□		□	□

• Standard □ Optional

Model		CM10	CM12	CM100	CM120	CM150
Software options						
Philips basic arrhythmia		•	•	•	•	•
Philips enhanced arrhythmia	C01	☐	☐	☐	☐	☐
Assisting venous puncture	C04	☐	☐	☐	☐	☐
Barcode reader support ⁶	C08	☐	☐	☐	☐	☐
Philips ST Map	C10		☐		•	•
48 hours of full disclosure (all waveforms)	C15	☐	☐	☐	☐	☐
HL7 output	C20	☐	☐	•	•	•
ADT inbound ⁷	C21			☐	☐	☐
Features and hardware options						
10-inch display		•		•		
12-inch display			•		•	
15-inch display						•
Touchscreen	E01	☐	☐	•	•	•
Maximum number of waveforms displayed		8	12	8	8	12
Bed-to-Bed overview		•	•	•	•	•
Early warning scoring (EWS)	C22	☐	☐	☐	☐	☐
External display connection		•	•	•	•	•
Internal recorder ⁸	E10	☐	☐	☐	☐	☐
Calculators (hemodynamics, Oxygenation, Drug, Renal)	C17	☐	☐	☐	☐	☐
Display mini trend, 8 hours			•		•	•
USB ports		2	2	2	2	2
Night mode		•	•	•	•	•
ECG analog output	E03	☐	•	•	•	•
Dual speaker				•	•	•
HL7 wireless LAN connectivity	E20	☐	☐	☐	☐	☐
Dashboard connectivity		•	•	•	•	•
Bed rail hook	E16	☐	☐	☐	☐	
Accessories and cable options						
Accessories not included	A02	☐	☐	☐	☐	☐
Adult accessories kit (3-lead) ^{9,10}	A01	☐	☐	☐	☐	☐
Pediatric accessories kit (3-lead) ^{10,11}	A03	☐	☐	☐	☐	☐
Neonatal accessories kits (3-lead) ^{10,12}	A04	☐	☐	☐	☐	☐
Adult premium accessories kit (3-lead) ^{2,10}	A05	☐	☐	☐	☐	☐
Adult premium accessories kit (5-lead) ^{2,10}	A06	☐	☐	☐	☐	☐
CO ₂ Mainstream sensor ¹³	A21	☐	☐	☐	☐	☐
Reusable adult/pediatric CO ₂ Mainstream adapter ¹⁴	K21	☐	☐	☐	☐	☐
Battery options						
3-cell lithium-ion battery	E11	☐	☐			
9-cell lithium-ion battery	E12	☐	☐	•	☐	☐
Dual 3-cell lithium-ion batteries	E13		☐ ¹⁵			
Dual 9-cell lithium-ion batteries	E14		☐ ¹⁵		☐	☐

• Standard ☐ Optional

The following are trademarks or registered trademarks of Masimo Corporation: DCI, FastSat, LNCS, LNOP, Masimo, PVI, SET, SpCO, SpHb, SpMet, TF-I, rainbow ReSposable, APOD, rainbow, M-LNCS, rainbow Acoustic Monitoring, rainbow Acoustic Sensor, rainbow Acoustic.

Oridion is a trademark and Microstream is a registered trademark of Oridion Medical Ltd.

Nellcor and OxiMax are trademarks of Medtronic.

1. Default unless Masimo or Nellcor selected.
2. May not be available in all countries.
3. Can only be ordered in combination with B02.
4. Only one of B20, B21 or B22 can be ordered.
5. Can only be ordered in combination with B10.
6. Barcode reader must be ordered separately:
- 98980 317 6611
- 98980 314 7821
7. Philips IntelliBridge Enterprise Interoperability solution required.
8. Can only be ordered together with either option E12 or option E14.
9. 3-lead ECG cable (grabber), adult SpO₂ sensor, adult NBP cuff and NBP hose.

10. SpO₂ sensor and cable only included when B01 is selected.
11. 3-lead ECG cable (grabber), pediatric SpO₂ sensor, pediatric NBP cuff and NBP hose.
12. 3-lead trunk cable, disposable neonatal ECG lead wires/electrodes, Philips neonatal hand/foot SpO₂ sensor, SpO₂ extension cable, neonatal NBP hose, 25 (5 cuffs per each size) neonatal single-patient NBP cuffs (various sizes).
13. Can only be ordered in combination with B21 Respironics Capnostat Mainstream sensor M2501A.
14. Can only be ordered in combination with A21.
15. Can only be ordered in combination with E01.

© 2021 Koninklijke Philips N.V. All rights are reserved. Philips reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.



philips.com

Printed in The Netherlands.
4522 991 70261 * JUL 2021

Philips Medizin Systeme
Böblingen GmbH
Hewlett-Packard Strasse 2
71034 Böblingen
Germany

