

PHILIPS

sonicare

Two sonic electric
toothbrushes

FlexCare

3 modes

4 brush heads

UV brush head sanitizer



HX6934/34



Advanced Cleaning

Innovative technology delivers better results. Now there's a guaranteed way to better oral health. It's the new Philips Sonicare FlexCare that adapts to your oral care needs.

Guides you to follow dental recommendations

- Two-minute timer helps ensure recommended brushing time
- Quadpacer interval timer encourages thorough brushing

Provides a superior clean

- Sonicare dynamic cleaning action drives fluid between teeth
- Angled brush head neck for better reach to back teeth
- UV brush head Sanitizer - Cleaner brush, better oral health

Proven to improve oral health

- Clinically proven safe and gentle
- Naturally whiter teeth with patented sonic technology
- Improves gum health in only two weeks

Offers a customized brushing experience

- Easy-start program to acclimate to the Sonicare experience
- Three different modes for a better cleaning experience



Highlights

Improves gum health



This Philips Sonicare electric toothbrush provides optimal cleaning between teeth and along the gumline for improved gum health in just two weeks.

Safe and gentle



Philips Sonicare is a gentle electric toothbrush for braces (brush heads wear out sooner when used on braces), and is safe for dental restorations (fillings, crowns, veneers) and periodontal pockets, too.

Sonic technology




Philips Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gumline.

SmarTimer




2-minute timer on this Philips Sonicare electric toothbrush helps ensure dental professional recommended brushing time

Angled brush head neck




The uniquely angled brush head neck of this toothbrush head makes it easier to reach back teeth, removing plaque in those hard-to-reach areas.

Quadpacer




30 seconds interval timer indicates when you have completed each quadrant of your mouth and signals for you to move on resulting in more consistent cleaning throughout your mouth

Naturally whiter teeth



The dynamic cleaning action of this toothbrush and the direct broader contact with each tooth, has proven to remove everyday stains, to help keep teeth naturally whiter.

Three different cleaning modes




The three modes include a clean mode for an outstanding clean in 2 minutes. A sensitive mode for gentle yet thorough cleaning of the gums and a massage mode to stimulate the gum.

UV brush head Sanitizer

All-in-one sanitizing, charging and storage. The UV brush head Sanitizer helps keep your brush head free of bacteria. It kills up to 99% of bacteria and viruses.*

Easy-start program



Gently increases power over first 14 uses to ease into the Philips Sonicare electric toothbrush experience

Specifications

Modes

Clean: For exceptional everyday clean

Sensitive: Gentle teeth and gum cleaning

Massage: Invigorating massage

2 Routines: Go Care and Max Care

Items included

Handles: 2 FlexCare

Brush heads: 2 InterCare standard, 1

ProResults standard, 1 ProResults compact

Charger: 1

UV sanitizer

Travel case: 2

Design and finishing

Color: White & green

Cleaning performance

Speed: Up to 62000 brush movement/min

Performance: Removes up to 83% of plaque in hard-to-reach areas

Health benefits: Helps improve gum health

Whitening benefits: Whitens teeth up to 2 shades

Timer: Quadpacer and SmarTimer

Ease of use

Brush head system: Easy click-on brush heads

Brushing time: Up to 2 weeks

Battery indicator: Light shows battery status

Handle: Slim ergonomic design, Rubber grip for easy handling

Technical specifications

Battery: Rechargeable

Battery type: Lithium ION

Operating time (full to empty): Up to 2 weeks

Power

Voltage: 110-220 V

Service

Warranty: 2-year limited warranty

